

Perkin Elmer - ICP 5000, 5500, 5500B, 6000, 6500 Glassware



- Provide an exact match for OEM Glassware
- Offer a complete product line
 - Save shipping cost
- Available in-stock and ready-to-ship



With every purchase of glassware, receive a **FREE PlasmaCAL ICP-AES/MS, 125 ml, 1000 µg/ml standard***

Torches & Accessories

Description	Cross Ref.	Catalog #
Torch, One Piece, w/Ball Joint	O047-2048	020-051-003
Torch Extension, Slotted	O047-3560	020-050-061
Torch Tube, Type II	N058-0528	020-050-063
Torch Tube, Type I, Demountable	N058-0233	020-050-064
Quartz Injector Tip, 1.2 mm for Type I Torch	---	020-051-301
Quartz Injector Tip, 1.6 mm for Type I Torch	N058-0449	020-050-067
Quartz Injector Tip, 2.0 mm for Type I Torch	---	020-050-010
Quartz Injector Tip, 1.2 mm for Type II Torch	N068-1631	020-050-182
Quartz Injector Tip, 1.6 mm for Type II Torch	O047-3292	020-050-183
Torch Extension	O047-3052	020-051-303
Torch Extension with Optical Port	---	020-051-305
Torch Extension, Slotted	O047-3560	020-051-307
Demountable Torch Bonnet, Type 1	O047-2033	020-051-309
Purge Tube	O047-3187	020-050-218
Purge Tube for Slotted Extension	O047-3561	020-050-217

Spray Chambers & Accessories

Description	Cross Ref.	Catalog #
Scott Spray Chamber, Glass	O047-2047	020-050-007
Cyclonic Spray Chamber	---	020-051-119
Cyclonic, Quartz	---	020-051-120
Cyclonic, Water Cooled	---	020-051-121
Cyclonic Spray, MLN	---	020-051-122
Cyclonic, 12/5 Socket, 4 mm Drain	---	020-060-124
Cyclonic, Water Cooled, 12/5 Socket, 4 mm Drain	---	020-051-125

*Contact SCP SCIENCE for more details.

Providing Innovative Solutions

to Analytical Chemists

Nebulizers

Description	Cross Ref.	Catalog #
Noordermeer	---	020-060-602
General (A3), concentric, 25-40 psi, 2.5-3.5ml/min	0047-2020	020-060-002
Hi-Solids (C1), concentric, 25-40 psi, 0.7-1.4ml/min	0047-2022	020-060-003
Hi-Solids (CO.5), concentric, 41-60 psi, 0.30.6ml/min	N810-2012	020-060-004
Hi-Solids (C3), concentric, 25-40 psi, 2.5-3.5ml/min	N810-2011	020-060-007
Low Flow (K3), concentric, 25-40 psi, 2.3-4.0ml/min	N068-1574	020-060-009
Teflon Mira Mist, concentric, 35-45 psi, 0.2-2.5ml/min	---	020-060-306
Peek Mira Mist, concentric, 45 psi, 0.2-2.5ml/min	N077-5330	020-060-308


Peristaltic Pump Tubing*

Description	Cross Ref.	Catalog #
Orange/Red, 0.0075", PVC, 2-stops, 12/pak	N069-5476	020-030-002
Black/Black, 0.030", PVC, 2-stops, 12/pak	0990-8587	020-030-007
Red/Red, 0.045", PVC, 2-stops, 12/pak	N068-0375	020-030-010
Black/White, 0.125", PVC, 2-stops, 12/pak	N812-2012	020-030-021
Black/Black, 0.030", PVC, 3-stops, 12/pak	B050-6058	020-033-007
White/White, 0.040", PVC, 3-stops, 12/pak	B050-6742	020-033-009
Red/Red, 0.045", PVC, 3-stops, 12/pak	B019-3160	020-033-010
Purple/Purple, 0.080", PVC, 3-stops, 12/pak	B019-9034	020-033-015
Yellow/Blue/Yellow, 0.060", PVC, 3-stops, 12/pak	B019-3161	020-033-019
Black/White/Black, 0.125", PVC, 3-stops, 12/pak	B050-8310	020-033-021
Black/Black, 0.030", Solvent Flex, 2-stops, 12/pak	0047-3550	020-030-107
Red/Red, 0.045", Solvent Flex, 2-stops, 12/pak	N068-0376	020-030-110
White/White, 0.040", Solvent Flex, 3-stops, 12/pak	B050-7692	020-033-109
Black/Black, 0.030", Silicon, 2-stops, 12/pak	N047-3552	020-030-207
Red/Red, 0.045", Silicon, 2-stops, 12/pak	N069-1595	020-030-210
Black/Black, 0.030", Viton, 2-stops, 12/pak	0990-8587	020-030-307
Perkin Elmer Mixing Kit, PVC, 3-stops, Orange/Orange:Orange/Green, 4/pak	---	020-049-602
Perkin Elmer Sciex Mixing Kit, PVC, 2-stops, Black/Black:Orange/Red, 4/pak	---	020-049-603

* For a complete list, please contact **SCP SCIENCE**

For more information or to place an order, please contact your
SCP SCIENCE Sales Representative or distributor.

USA
Tel.: (800) 361-6820
Fax: (800) 253-5549



Canada / International
Tel.: (800) 361-6820 / (514) 457-0701
Fax: (800) 253-5549 / (514) 457-4499



Europe
Tel.: +33 (0)1 69 18 71 17
Fax: +33 (0)1 60 92 05 67