

Thermo-TJA Glassware

- Provide an exact match for OEM Glassware
- Offer a complete product line
 - Save shipping cost
- Available in-stock and ready-to-ship



With every purchase of glassware, receive a **FREE PlasmaCAL ICP-AES/MS, 125 ml, 1000 µg/ml standard***

Torches and Accessories

Description

Cross Ref.

Catalog #

2000,
4000

Polyscan,
Atomscan, Radial,
61e, Iris, POEMS

975, 9000,
1100

Trace Advantage
Axial, TraceScan,
Iris Ap

Description	Cross Ref.	Catalog #	2000, 4000	Polyscan, Atomscan, Radial, 61e, Iris, POEMS	975, 9000, 1100	Trace Advantage Axial, TraceScan, Iris Ap
Torch	19040073	020-050-112	•			
Torch bonnet (2)	19040075	020-049-300	•		•	
Ceramic base torch, high/ low flow w/ Injector 1.5 mm	126432-02	020-050-123		•		
Ceramic base torch, high flow w/ Injector 1.5 mm	126432-01	020-050-124		•		
Ceramic base torch, high flow w/ Injector 1.0 mm	125124-01	020-050-125		•		
Ceramic base torch, high flow w/ Alumina Inj. 1.5mm	127179-01	020-049-000		•		
Argon connectors ceramic base torches	N/A	020-051-200		•		
Injector tip, quartz, 1.0 mm	125407-01	020-050-118		•		•
Injector tip, quartz, 1.5 mm	125407-00	020-050-119		•		•
Injector tip, quartz, 2.0 mm	136502-01	020-051-340		•		•
Injector tip, quartz, 2.4 mm	136502-02	020-051-344		•		•
Injector tip, alumina, 1.0 mm	N/A	020-049-302		•		•
Injector tip, alumina, 1.5 mm	136436-01	020-051-300		•		•
Anti-Arc cap	085207-00	020-050-126		•		
TJA anti-arc cap/ Flat edge	085207-01	020-050-127		•		
Torch w/ removable 1.0 mm top	19040054	020-050-110			•	
As above, but with 1 extra 1.0 mm inj.	N/A	020-050-111			•	
As above, but with 1 extra 1.25 mm inj.	N/A	020-050-195			•	
As above, but with 1 extra 1.5 mm inj.	N/A	020-050-196			•	
As above, but with 1 extra 1.7 mm inj.	N/A	020-049-196			•	
HF/CR Torch body for alum./sapphire Inj.	N/A	020-050-193			•	
Torch extension	N/A	020-050-117			•	
Torch extension w/ optical port	N/A	020-049-301			•	
Ceramic base torch w/ 1.5 mm Inj. tip	126432-03	020-050-192				•
Ceramic base torch - DV w/ 1.5 mm Inj. tip	138562-00	020-051-192				•
Optic sleeve / port DV	138563-00	020-049-304				•

*Contact **SCP SCIENCE** for more details.

Spray Chambers and Accessories

Description

Cross Ref.

Catalog

2000,
4000

Polyscan,
Atomscan, Radial,
61e, Iris, POEMS

975, 9000,
1100

Trace Advantage
Axial, TraceScan,
Iris AP

Standard	12794700	020-050-114	•			
Aerosol baffle for 020-050-114	00395200	020-050-069	•		•	
Cyclonic	N/A	020-050-124	•			
Teflon coated viton O-Ring	N/A	020-050-043	•			
Curved drain, open-ended	13008900	020-051-116		•		
Aerosol baffle for spray chamber	00395200	020-050-069		•		
Cyclonic, with 12 mm drain	13610500	020-051-108		•		
Cyclonic, water jacket, with 12 mm drain	N/A	020-051-109		•		
Standard	19040065	020-050-113			•	
Curve drain	N/A	020-050-115			•	
Curve drain, open ends, water cooled	N/A	020-050-190			•	
Cyclonic, 2 piece w/ 12 mm drain	13612000	020-051-112				
Cyclonic, w/ 12 mm drain	N/A	020-050-191				•
Cyclonic, water jacket, w/ 12 mm drain	N/A	020-049-101				•

Nebulizers

Description

Cross Ref.

Catalog

2000,
4000

Polyscan,
Atomscan, Radial,
61e, Iris, POEMS

975, 9000,
1100

Trace Advantage
Axial, TraceScan,
Iris AP

Cross flow	126989-00	020-060-606	•	•	•	
Noordermeer	N/A	020-060-602		•		•
General (A2), concentric, 25-40 psi, 1.5-2.4ml/min	N/A	020-060-001		•		•
Hi-Solids (C2), concentric, 25-40 psi, 1.5-2.4ml/min	N/A	020-060-008		•		•
Low Flow (K2), concentric, 25-40 psi, 1.0-2.2ml/min	N/A	020-060-006		•		•
Teflon Mira Mist, concentric, 35-45 psi, 0.2-2.5ml/min	136388-00	020-060-306		•		•
Peek Mira Mist, concentric, 45 psi, 0.2-2.5ml/min	N/A	020-060-308		•		•

3-Stop Peristaltic Pump Tubing*

Description

Cross Ref.

Catalog

2000,
4000

Polyscan,
Atomscan, Radial,
61e, Iris, POEMS

975, 9000,
1100

Trace Advantage
Axial, TraceScan,
Iris AP


Orange/Red/Orange, 0.0075", PVC, 12/pak	N/A	020-033-002	•	•	•	•
Orange/Green/Orange, 0.015", PVC, 12/pak	130177-06	020-033-004	•	•	•	•
Red/Red, 0.045", PVC, 12/pak	130117-01	020-033-010	•	•	•	•
Orange/Orange, 0.035", Silicon, 6/pak	130177-02	020-030-008	•	•	•	•
Red/Red, 0.045", Silicon, 6/pak	130117-03	020-033-208	•	•	•	•
Orange/Orange/, 0.035", Viton, 12/pak	130177-04	020-033-210	•	•	•	•
Red/Red, 0.045", Viton, 12/pak	130177-05	020-033-328	•	•	•	•
Internal Standard Mixing Kit, 4/pak	136708-00	020-049-601	•	•	•	•

*For a complete list, please contact **SCP SCIENCE**.

For more information or to place an order, please contact your
SCP SCIENCE Sales Representative or distributor.


USA

Tel.: (800) 361-6820
Fax: (800) 253-5549



Canada / International

Tel.: (800) 361-6820 / (514) 457-0701
Fax: (800) 253-5549 / (514) 457-4499



Europe

Tel.: +33 (0)1 69 18 71 17
Fax: +33 (0)1 60 92 05 67