

Varian Liberty Series and Vista Glassware

- Provide an exact match for OEM Glassware
- Offer a complete product line
 - Save shipping cost
- Available in-stock and ready-to-ship



Torches & Accessories - Radial



Description	Cross Ref.	Catalog #
Torch, Low Flow, ICP-AES	20-100696-00 / 20-100696-90	020-050-140
Torch, Varian ICP-AES with Ball Joint	20-100696-00B	020-050-002
Torch, Extended for Nitrogen Determination	20-100764-00	020-050-001
Torch Kit, ICP-AES Demountable, HF resistant	99-100568-00	020-050-130
Demountable Torch Body Only, HF resistant	20-100729-90	020-050-142
Alumina Injector - 1.8 mm	24-100675-90	020-050-145
Alumina Injector - 2.5 mm	---	020-052-346
Quartz Injector - 1.8 mm	99-100570-00 (3/pack) / 20-100730-00	020-050-144
Quartz Injector - 0.8 mm	20-100775-00 (3/pack)	020-052-301
Teflon Injector Holder, Varian	99-100573-00	020-050-199
Torch Bonnet (pack of 2)	20-100707-90	020-050-143

Torches & Accessories - Axial



Description	Cross Ref.	Catalog #
Torch, Low Flow, One Piece	20-100812-00	020-052-002
Torch, Low Flow, One Piece w/90 Degree Stem Bend	20-100904-00	020-054-125
Torch, Low Flow, One Piece, High Solids	20-100948-00	020-054-080
Torch Kit, Demountable, HF Resistant	99-100847-00	020-052-003
Demountable Torch Body Only	20-100834-00	020-052-300
Alumina Injector - 1.8 mm	24-100675-90	020-050-145
Alumina Injector - 2.5 mm	---	020-052-346
Quartz Injector - 1.8 mm	99-100570-00 (3/pack) / 20-100730-00	020-050-144
Quartz Injector - 0.8 mm	20-100775-00 (3/pack)	020-052-301
Varian Teflon Injector Holder	99-100573-00	020-050-199

*Contact **SCP SCIENCE** for more details.

Spray Chambers & Accessories



Description	Cross Ref.	Catalog #
Cyclonic Spray Chamber, 12/5 Socket Top, 4 mm Drain	---	020-060-124
Optimization solution, Mn, 5 ppm, 125 ml	---	140-130-251

Nebulizers



Description	Cross Ref.	Catalog #
Concentric (C-type, 1 ml/min Ar @ 35 psi)	20-100765-00 / 20-100816-00	020-060-003
Concentric (C-type, 2.5-3.5 ml/min Ar @ 25-40 psi), High Solids	20-100765-00	020-060-006
Concentric (K-type, 1.0-2.2 ml/min Ar @ 25-40 psi), Low Flow	20-100695-00	020-060-008
Concentric (K-type, 2.3-4.0 ml/min Ar @ 25-40 psi), Low Flow	20-100816-00	020-060-009
Noordermeer V-Groove MDSN Nebulizer, 0.7 L/min to 1.0 @ 35 psi	---	020-060-602
High Solids V-Groove Nebulizer, 0.7 lpm Ar @ 35 psi	99-100574-00	020-060-607
Burgener MiraMist, 0.2 - 2.5 ml/min @ 45 psi (Teflon)	---	020-060-306
Burgener PEEK MiraMist, 0.2 - 2.5 ml/min @ 45 psi (PEEK)	---	020-060-308


Peristaltic Pump Tubing*

Description	Cross Ref.	Catalog #
Black/Black, 0.030", PVC, 2-stops, 12/pak	37-100272-00	020-030-007
White/White, 0.040", PVC, 2-stops, 12/pak	37-100344-00	020-030-009
Grey/Grey, 0.050", PVC, 2-stops, 12/pak	37-100345-00	020-030-011
Blue/Blue, 0.065", PVC, 2-stops, 12/pak	37-100346-00	020-030-013
Purple/Black, 0.090", PVC, 2-stops, 12/pak	37-100273-00	020-030-016
Black/Black, 0.030", Solvent Flex, 2-stops, 12/pak	37-100348-00	020-030-107
White/White, 0.040", Solvent Flex, 2-stops, 12/pak	37-100350-00	020-030-109
Grey/Grey, 0.050", Solvent Flex, 2-stops, 12/pak	37-100352-00	020-030-111


*For a complete list and to inquire about your FREE ICP-AES/MS Standard, please contact **SCP SCIENCE**.

For more information or to place an order, please contact your
SCP SCIENCE Sales Representative or distributor.

USA
Tel.: (800) 361-6820
Fax: (800) 253-5549



Canada / International
Tel.: (800) 361-6820 / (514) 457-0701
Fax: (800) 253-5549 / (514) 457-4499



Europe
Tel.: +33 (0)1 69 18 71 17
Fax: +33 (0)1 60 92 05 67