


Agilent -

ICP-MS Supplies for HP 4500

- Provide an exact match for Agilent 4500
- Offer a complete product line
 - Series of Nebulizers, Torches, & Spray Chambers
- Ensure accuracy and peace of mind



With every purchase of glassware, receive a **FREE PlasmaCAL** ICP-AES/MS, 125 ml, 1000 µg/ml standards*

Torches & Accessories



Description	Cross Ref.	Catalog #
Torch with 2.5 mm Injector	G1820-65028	020-050-139
Torch, HF Demountable	G1820-65035	020-050-135
HP Teflon Injector Holder	G1820-65015	020-050-136

Spray Chambers & Accessories




Description	Cross Ref.	Catalog #
Quartz spray chamber for Babington Nebulizer	G1820-65337	020-056-101
Quartz spray chamber for Cross Flow Nebulizer	G1820-65339	020-056-102
Quartz ball joint connector	G1820-80237	020-056-201
Clamps	G1820-65007	020-048-201

Nebulizers



Description	Cross Ref.	Catalog #
Babington Nebulizer	G1820-60453	020-060-605
Pyrex Concentric Nebulizer	G1820-65030	020-060-002
Quartz Concentric Nebulizer	G1820-65138	020-060-011
Cross Flow Nebulizer	G1833-65462	020-060-604

*Contact **SCP SCIENCE** for more details.



SCP SCIENCE provides the exact matched Cones for Agilent 4500. *PlasmaFORM* ICP-MS Cones are manufactured by the electroforming process rather than by traditional CNC lathe methods. Providing better performance with:


- Sharper hole edges and superior surface finishes
- Minimize gas disturbances
- Enhanced thermal characteristics
- Extended lifetimes

ICP-MS Cones

Description	Cross Ref.	Catalog #
Nickel Sampling	G1820-65238	020-107-051
Nickel Skimmer	G1820-65265	020-107-052
Nickel Sampling, for T-Mode	G1820-65480	020-107-055
Nickel Skimmer, for T-Mode	G1820-65481	020-107-056
Platinum Sampling, copper base	G1820-65239	020-107-151
Platinum Skimmer, copper base	G1820-65237	020-107-152

Peristaltic Pump Tubing

Tubing Type	Color Code	Stops	ID mm	Quantity	Catalog #
PVC	Orange/Blue	2	0.25	12	020-030-003
PVC	White/White	2	1.02	12	020-030-009
PVC	Orange/Blue	3	0.25	12	020-033-003
PVC	Orange/Orange	3	0.89	12	020-033-008
Solvent Flex	Orange/Blue	2	0.25	12	020-030-103
Solvent Flex	Orange/Orange	2	0.89	12	020-030-108
Solvent Flex	White/White	2	1.02	12	020-030-109
Solvent Flex	Orange/Blue	3	0.25	12	020-033-103
Solvent Flex	Orange/Orange	3	0.89	12	020-033-108
Solvent Flex	White/White	3	1.02	12	020-033-109
Silicon	Orange/Orange	2	0.89	6	020-030-208
Silicon	White/White	2	1.02	6	020-030-209
Silicon	Orange/Orange	3	0.89	6	020-033-208
Silicon	White/White	3	1.02	6	020-033-209
Viton	Orange/Orange	2	0.89	12	020-030-308
Viton	White/White	2	1.02	12	020-030-309
Viton	Orange/Orange	3	0.89	12	020-033-308
Viton	White/White	3	1.02	12	020-033-309



*For a complete list, please contact **SCP SCIENCE**.

For more information or to place an order, please contact your
SCP SCIENCE Sales Representative or distributor.

USA

Tel.: (800) 361-6820
Fax: (800) 253-5549

MK-MKG003-45AG-1.0-E

Canada / International

Tel.: (800) 361-6820 / (514) 457-0701
Fax: (800) 253-5549 / (514) 457-4499

Europe

Tel.: +33 (0)1 69 18 71 17
Fax: +33 (0)1 60 92 05 67