

—Varian ICP-MS & Ultramass Glassware—

SCP SCIENCE maintains a full range of glassware to keep your analyses running.

- Provide an exact match for OEM Glassware
- Offer a complete product line
 - Save shipping cost
- Available in-stock and ready-to-ship



With every purchase of nebulizer, spray chamber, torch or cone, receive a **FREE PlasmaCAL ICP-AES/MS, 125 ml, 1000 µg/ml standard***

Torches & Accessories



020-050-198

Description	Cross Ref.	Catalog #
Torch, One Piece, Low Flow	20-100783-00	020-050-141
Torch, Demountable, Body Only	20-100793-00	020-050-198
Torch, One Piece	20-101007-00	020-054-308
Demountable Torch, Complete	99-100717-00	020-052-004
Quartz Injector, 0.8mm	20-100775-00	020-052-301
Quartz Injector, 1.8mm	20-100730-00	020-050-144
Alumina Injector 1.8mm	24-100675-00	020-050-145
Alumina Injector 2.5mm	---	020-052-346
Teflon Injector Holder	99-100573-00	020-050-199
Sheath Gas, Bypass	20-101008-00	020-054-304
Sheath Gas, Port	20-101009-00	020-054-306
Torch Bonnet	20-100707-00	020-050-143

Spray Chambers



020-050-174

Description	Cross Ref.	Catalog #
Cyclonic Spray Chamber, 12/5 Socket Top, 4mm Drain	---	020-050-174
Cyclonic Spray Chamber	20-101010-00	020-054-310

*Contact SCP SCIENCE for more details.

Nebulizers




Description	Cross Ref.	Catalog #
Concentric (C-type, 1 ml/min Ar @ 35 psi)	20-100816-00	020-060-003
Concentric (C-type, 2.5-3.5 ml/min Ar @ 25-40 psi), High Solids	20-100765-00	020-060-006
Concentric (K-type, 1.0-2.2 ml/min Ar @ 25-40 psi), Low Flow	20-100695-00	020-060-008
Concentric (K-type, 2.3-4.0 ml/min Ar @ 25-40 psi), Low Flow	20-100816-00	020-060-009
High Solids V-Groove Nebulizer, 0.7 lpm Ar @ 35 psi	99-100574-00	020-060-607
Babbington Glass V-Groove, 1.0 L/min @ 35 psi	---	020-060-605
Burgener MiraMist, 0.2 - 2.5 ml/min @ 45 psi (Teflon)	---	020-060-306
Burgener PEEK MiraMist, 0.2 - 2.5 ml/min @ 45 psi (PEEK)	---	020-060-308

ICP-MS Cones



Description	Cross Ref.	Catalog #
Nickel Sampler F-Series Type	99-100940-00	020-107-070
Nickel Skimmer F-Series Type	99-100941-00	020-107-072

Peristaltic Pump Tubing*



Description	Cross Ref.	Catalog #
Black/Black, 0.030", PVC, 2-stops, 12/pak	37-100272-00	020-030-007
White/White, 0.040", PVC, 2-stops, 12/pak	37-100344-00	020-030-009
Grey/Grey, 0.050", PVC, 2-stops, 12/pak	37-100345-00	020-030-011
Blue/Blue, 0.065", PVC, 2-stops, 12/pak	37-100346-00	020-030-013
Purple/Black, 0.090", PVC, 2-stops, 12/pak	37-100273-00	020-030-016
Black/Black, 0.030", Solvent Flex, 2-stops, 12/pak	37-100348-00	020-030-107
White/White, 0.040", Solvent Flex, 2-stops, 12/pak	37-100350-00	020-030-109
Grey/Grey, 0.050", Solvent Flex, 2-stops, 12/pak	37-100352-00	020-030-111
Black/Black, 0.030", Silicon, 2-stops, 6/pak	37-100441-00	020-030-207
White/White, 0.040", Silicon, 2-stops, 6/pak	37-100442-00	020-030-209
Grey/Grey, 0.050", Silicon, 2-stops, 6/pak	37-100443-00	020-030-211
Blue/Blue, 0.065", Silicon, 2-stops, 6/pak	37-100444-00	020-030-213
Black/Black, 0.030", Viton, 2-stops, 12/pak	37-100272-00	020-030-307
White/White, 0.040", Viton, 2-stops, 12/pak	37-100344-00	020-030-309
Grey/Grey, 0.050", Viton, 2-stops, 12/pak	37-100345-00	020-030-311
Blue/Blue, 0.065", Viton, 2-stops, 12/pak	37-100346-00	020-030-313

*For a complete list and to inquire about your FREE ICP-AES/MS Standard, please contact **SCP SCIENCE**.

For more information or to place an order, please contact your
SCP SCIENCE Sales Representative or distributor.

USA
Tel.: (800) 361-6820
Fax: (800) 253-5549

Canada / International
Tel.: (800) 361-6820 / (514) 457-0701
Fax: (800) 253-5549 / (514) 457-4499

Europe
Tel.: +33 (0)1 69 18 71 17
Fax: +33 (0)1 60 92 05 67