

# DigiFILTER

1, 2, 3...

**PARTICULATE FREE!**

*The Next Generation in filtering samples has arrived.*

*The vacuum assisted DigiFILTER assembly provides a quick and easy way to filter samples prior to analysis.*



*\*patent pending*

**SCP SCIENCE**

[www.scpscience.com](http://www.scpscience.com)

**USA**

Tel.: (800) 361-6820  
Fax: (800) 253-5549

**Canada / International**

Tel.: (800) 361-6820 / (514) 457-0701  
Fax: (800) 253-5549 / (514) 457-4499

**Europe**


Tel.: +33 (0)1 69 18 71 17  
Fax: +33 (0)1 60 92 05 67

## 4, 5, 6... Disposable, easy and quick

# DigiFILTER

Filter 12 samplers at a time using the DigiFILTER manifold

Our patent-pending **DigiFILTERs** will increase lab productivity by filtering particulates from digested samples in a matter of seconds. Available in 0.45 or 1.0  $\mu\text{m}$  hydrophilic Teflon® membranes.



Item	Catalog Number
DigiFILTER, 0.45 micron Teflon® (pk/100)	010-500-070
DigiFILTER, 1.0 micron Teflon® (pk/100)	010-500-080

Item	Catalog Number
DigiFILTER Manifold	010-500-230
18 Position Rack	010-500-025

### DigiFILTER Accessories

Item	Cat. Number
DigiTUBEs, RackLock (750/Case)	010-500-261
DigiTUBEs, Non-RackLock (750/Case)	010-500-263
Vacuum Pump, Standard (115V)	010-500-235
Vacuum Pump, Teflon® lined (115V)	010-790-530
Vacuum Pump, Standard (230V)	010-790-509
Vacuum Pump, Teflon® lined (230V)	010-790-531



Catalogue Number:  
010-500-070

Catalogue Number:  
010-500-080