

Atomic Absorption Calibration Standards

SCP SCIENCE has a wide selection of Standards; more than has for Flame and Graphite Furnace Atomic Absorption

Satisfaction Guaranteed

"For best results, use SCP SCIENCE AA Single element Calibration Standards"



- Certificate of Analysis with actual matrix, concentration and traceability to NIST
- Complete documentation for audit purposes
- 2 expiry dates (up to 21 months unopened & 15 months opened)
- Immediate availability for most common elements

Catalog Number

Element	Symbol	Matrix	Code	1000 µg/ml 125 ml	1000 µg/ml 500 ml
Aluminum	Al	HCl	✓ ⑧	140-002-131	140-002-135
Antimony	Sb	HNO ₃ / tr. Tartaric Acid	✓ ⑧	140-001-511	140-001-515
Arsenic	As	HNO ₃	✓ ⑧	140-001-331	140-001-335
Barium	Ba	HNO ₃	✓ ⑧	140-001-561	140-001-565
Beryllium	Be	HNO ₃	✓ ⑧	140-001-041	140-001-045
Bismuth	Bi	HNO ₃	✓ ⑧	140-001-831	140-001-835
Boron	B	H ₂ O		140-000-051	140-000-055
Cadmium	Cd	HNO ₃	✓ ⑧	140-001-481	140-001-485
Calcium	Ca	HNO ₃	✓ ⑧	140-001-201	140-001-205
Chromium	Cr	HCl	✓ ⑧	140-002-241	140-002-245
Cobalt	Co	HNO ₃	✓ ⑧	140-001-271	140-001-275
Copper	Cu	HNO ₃	✓ ⑧	140-001-291	140-001-295
Gold	Au	HCl	✓ ⑧	140-002-791	140-002-795*
Iron	Fe	HNO ₃	✓ ⑧	140-001-261	140-001-265
Lead	Pb	HNO ₃	✓ ⑧	140-001-821	140-001-825
Lithium	Li	HNO ₃	✓ ⑧	140-001-031	140-001-035
Magnesium	Mg	HNO ₃	✓ ⑧	140-001-121	140-001-125
Manganese	Mn	HNO ₃	✓ ⑧	140-001-251	140-001-255
Mercury	Hg	HNO ₃	✓ ⑧	140-001-801	140-001-805
Molybdenum	Mo	H ₂ O		140-000-421	140-000-425
Nickel	Ni	HNO ₃	✓ ⑧	140-001-281	140-001-285
Potassium	K	HNO ₃	✓ ⑧	140-001-191	140-001-195
Selenium	Se	HNO ₃	✓ ⑧	140-001-341	140-001-345
Silicon	Si	H ₂ O / tr. HF		140-000-141	140-000-145
Silver	Ag	HNO ₃	✓ ⑧	140-001-471	140-001-475
Sodium	Na	HNO ₃	✓ ⑧	140-001-111	140-001-115
Strontium	Sr	HNO ₃	✓ ⑧	140-001-381	140-001-385
Tin	Sn	HCl	✓ ⑧	140-002-501	140-002-505
Titanium	Ti	H ₂ O / tr. HF		140-000-221	140-000-225
Vanadium	V	HNO ₃	✓ ⑧	140-001-231	140-001-235
Zinc	Zn	HNO ₃	✓ ⑧	140-001-301	140-001-305

Free Standards with Lamps & GFAA Tubes*

For a complete listing of SCP SCIENCE products, get your free Catalogue Today!!



SCP SCIENCE
www.scpscience.com

Sample Prep

ICP-AES & MS

AA

XRF

Calibration Standards

Certified Reference Materials

***Free Standards**



With every purchase of any Hollow Cathode Lamp or package of GFAA tubes, receive a 125 ml, 1000 µg/ml standard at no cost !!!

*Hazardous shipment charges are not included. The following element is excluded from all special offers: Au.

MK-MKG003-AAST-1.0-E

USA
Tel.: (800) 361-6820
Fax: (800) 253-5549

Canada / International
Tel.: (800) 361-6820 / (514) 457-0701
Fax: (800) 253-5549 / (514) 457-4499

Europe
Tel.: +33 (0)1 69 18 71 17
Fax: +33 (0)1 60 92 05 67