

SCP SCIENCE presents *Digi*TUBEs

Introducing the latest design for disposable digestion tubes;

The RackLock *Digi*TUBEs:

- * Allows one-handed closure of screw cap digestion tubes
- * Permits dilution to volume and tube closure in an easy two-handed operation
- * Saves sample preparation time



- * Accurate 50 ml graduation meeting Class A specifications for dilutions
- * Ultra low leachable metal content to avoid metal contamination
- * Leak-free screw cap to prevent spills
- * Used with *Digi*PREP systems, or with other digestion systems using RackLock Racks



New Low Volume Pricing

Save up to 16% when you purchase 5 or more cases (750 tubes/each) of *Digi*TUBEs (Cat.# 010-500-261 or 010-500-263).

Ordering Information


Catalog No.	Description	Quantity
010-500-261	<i>Digi</i> TUBE - Digestion Tube with RackLock	750
010-500-263	<i>Digi</i> TUBE - Digestion Tube without RackLock	750
010-500-025	18 Position RackLock Rack, for use with HotBlock™	1
010-500-021	24 Position RackLock Rack, for use with <i>Digi</i> PREP (48 positions)	1
010-505-021	24 position RackLock Rack, for use with <i>Digi</i> PREP Jr. (24 positions)	1
010-500-081	Disposable Polypropylene Watch Glasses	1000
010-500-028	<i>Digi</i> TUBE Foam Storage Racks	5

*Price per cases when 5 or more cases are purchased. HotBlock™ is a trademark of Environmental Express


USA
Tel.: (800) 361-6820
Fax: (800) 253-5549



Canada / International
Tel.: (800) 361-6820 / (514) 457-0701
Fax: (800) 253-5549 / (514) 457-4499



Europe
Tel.: +33 (0)1 69 18 71 17
Fax: +33 (0)1 60 92 05 67



Sample Prep
ICP-AES & MS
AA
XRF
Calibration Standards
Certified Reference Materials

Have you seen

*Digi*PREP HT
High Temperature
Digestion System




The *Digi*PREP HT is a new System for Kjeldahl digestion

- Versatile heating to 450°C
- Acid resistant graphite block
- Available in two sizes: 40 tube (100ml) or 20 tube (250ml) sample capacity

Call for more information or to request a quotation!!

More than a digestion tube...

The revolutionary *DigiTUBE* replaces:



- * Volumetric flasks
- * Graduated cylinders
- * Autosampler tubes
- * Storage containers
- * Digestion tubes

This single use, all-in-one universal vessel saves time, money, and eliminates the need to stock 5 different vessels!

Also Available for the Cleanroom

DigiTUBEs are available in shrinkwrap packaging (cases of 500) to ensure minimum contamination from storage and handling.

Ordering Information (continuation)

Catalogue No.	Description	Quantity
010-500-061	<i>DigiTUBE</i> - Digestion Tube with RackLock	500
010-500-062	<i>DigiTUBE</i> - Digestion Tube with RackLock, Shrink Wrap	20x25
010-500-063	<i>DigiTUBE</i> - Digestion Tube w/o RackLock	500
010-500-064	<i>DigiTUBE</i> - Digestion Tube w/o RackLock, Shrink Wrap	20x25



Did you know about the *DigiTUBE* Calibration Certificate?

Every case of *DigiTUBEs* comes complete with a Calibration Certificate that guarantees Class A tolerance at the 50 ml line of every *DigiTUBE*. Quality, consistency, and cost effectiveness... that's the **SCP SCIENCE** way.

SCP SCIENCE

www.scpscience.com

sales@scpscience.com