

— Perkin Elmer - Optima 2000 & 4000 Glassware —



- Provide an exact match for OEM Glassware
- Offer a complete product line
 - Save shipping cost
- Available in-stock and ready-to-ship



With every purchase of glassware, receive a **FREE** PlasmaCAL ICP-AES/MS, 125 ml, 1000 µg/ml standards*

Torches & Accessories

Description	Cross Ref.	Catalog #
Torch, DV (no injector)	N077-0338	020-051-002
Torch Body, Radial	N077-1500	020-051-005
Torch Body, Short Version	N077-0344	020-051-098
Quartz Injector, 2.0 mm	N077-5014	020-051-302
Quartz Injector, 3.0 mm	N077-5224	020-051-304
Quartz Injector, 1.6 mm	N077-5223	020-051-306
Quartz Injector, 1.2 mm	N077-5226	020-051-308
Quartz Injector, 0.8 mm	N077-5225	020-051-310
Alumina Injector, 0.8 mm	N077-5227	020-051-317
Alumina Injector, 1.2 mm	N077-5228	020-051-318
Alumina Injector, 2.0 mm	N077-5177	020-051-319
Torch Bonnet (pk/2)	N077-5289	020-051-350
Radial Purge Window, 4000	N077-0322	020-051-256
Axial Purge Window, 4000	N077-1116	030-050-158
Purge Tube Window, Radial, Short	N069-0672	020-050-254
Copper Strips (10)	N077-5297	020-051-001

Spray Chambers & Accessories

Description	Cross Ref.	Catalog #
Cyclonic, No Baffle	N077-6052	020-051-114
Cyclonic, Baffle	N077-6053	020-051-113
Mini Cyclonic	---	020-051-107
Cyclonic, 4 mm Baffle	N077-6090	020-051-138

*Contact SCP SCIENCE for more details.

Providing Innovative Solutions

to Analytical Chemists

Nebulizers

Description	Cross Ref.	Catalog #
Noordermeer	---	020-060-602
General (A3), concentric, 25-40 psi, 2.5-3.5ml/min	0047-2020	020-060-002
Hi-Solids (C1), concentric, 25-40 psi, 0.7-1.4ml/min	0047-2022	020-060-003
Hi-Solids (CO.5), concentric, 41-60 psi, 0.30.6ml/min	N810-2012	020-060-004
Hi-Solids (C3), concentric, 25-40 psi, 2.5-3.5ml/min	N810-2011	020-060-007
Low Flow (K3), concentric, 25-40 psi, 2.3-4.0ml/min	N068-1574	020-060-009
Teflon Mira Mist, concentric, 35-45 psi, 0.2-2.5ml/min	---	020-060-306
Peek Mira Mist, concentric, 45 psi, 0.2-2.5ml/min	N077-5330	020-060-308


Peristaltic Pump Tubing*

Description	Cross Ref.	Catalog #
Orange/Red, 0.0075", PVC, 2-stops, 12/pak	N069-5476	020-030-002
Black/Black, 0.030", PVC, 2-stops, 12/pak	0990-8587	020-030-007
Red/Red, 0.045", PVC, 2-stops, 12/pak	N068-0375	020-030-010
Black/White, 0.125", PVC, 2-stops, 12/pak	N812-2012	020-030-021
Black/Black, 0.030", PVC, 3-stops, 12/pak	B050-6058	020-033-007
White/White, 0.040", PVC, 3-stops, 12/pak	B050-6742	020-033-009
Red/Red, 0.045", PVC, 3-stops, 12/pak	B019-3160	020-033-010
Purple/Purple, 0.080", PVC, 3-stops, 12/pak	B019-9034	020-033-015
Yellow/Blue/Yellow, 0.060", PVC, 3-stops, 12/pak	B019-3161	020-033-019
Black/White/Black, 0.125", PVC, 3-stops, 12/pak	B050-8310	020-033-021
Black/Black, 0.030", Solvent Flex, 2-stops, 12/pak	0047-3550	020-030-107
Red/Red, 0.045", Solvent Flex, 2-stops, 12/pak	N068-0376	020-030-110
White/White, 0.040", Solvent Flex, 3-stops, 12/pak	B050-7692	020-033-109
Black/Black, 0.030", Silicon, 2-stops, 12/pak	N047-3552	020-030-207
Red/Red, 0.045", Silicon, 2-stops, 12/pak	N069-1595	020-030-210
Black/Black, 0.030", Viton, 2-stops, 12/pak	0990-8587	020-030-307
Perkin Elmer Mixing Kit, PVC, 3-stops, Orange/Orange:Orange/Green, 4/pak	---	020-049-602
Perkin Elmer Sciex Mixing Kit, PVC, 2-stops, Black/Black:Orange/Red, 4/pak	---	020-049-603

* For a complete list, please contact **SCP SCIENCE**

For more information or to place an order, please contact your
SCP SCIENCE Sales Representative or distributor.

USA
Tel.: (800) 361-6820
Fax: (800) 253-5549



Canada / International
Tel.: (800) 361-6820 / (514) 457-0701
Fax: (800) 253-5549 / (514) 457-4499



Europe
Tel.: +33 (0)1 69 18 71 17
Fax: +33 (0)1 60 92 05 67