

— Perkin Elmer - Optima 3000 & 3000SC Glassware —



- Provide an exact match for OEM Glassware
- Offer a complete product line
 - Save shipping cost
- Available in-stock and ready-to-ship



With every purchase of glassware, receive a **FREE** PlasmaCAL ICP-AES/MS, 125 ml, 1000 µg/ml standards*

Torches & Accessories

Description	Cross Ref.	Catalog #
Torch Body	N069-0568	020-050-181
Quartz Sample Injector 1.2 mm	N068-1631	020-050-182
Quartz Sample Injector 1.6 mm	N069-5448	020-050-238
Quartz Sample Injector 2.0 mm	N069-5441	020-050-234
Torch Bonnet, Flat Edge (pk/2) (before Sept. 1994)	N069-5295	020-050-185
Torch Bonnet (pk/2) (after Sept. 1994)	N069-5456	020-050-204
Purge Extension Window	N058-1497	020-050-213

Spray Chambers & Accessories

Description	Cross Ref.	Catalog #
Scott, Type II	---	020-051-115
Cyclonic, Optima	---	020-051-116
Cyclonic, for Modified Lichte Nebulizer	---	020-060-126

*Contact SCP SCIENCE for more details.

Providing Innovative Solutions

to Analytical Chemists

Nebulizers

Description	Cross Ref.	Catalog #
Noordermeer	---	020-060-602
General (A3), concentric, 25-40 psi, 2.5-3.5ml/min	0047-2020	020-060-002
Hi-Solids (C1), concentric, 25-40 psi, 0.7-1.4ml/min	0047-2022	020-060-003
Hi-Solids (CO.5), concentric, 41-60 psi, 0.30.6ml/min	N810-2012	020-060-004
Hi-Solids (C3), concentric, 25-40 psi, 2.5-3.5ml/min	N810-2011	020-060-007
Low Flow (K3), concentric, 25-40 psi, 2.3-4.0ml/min	N068-1574	020-060-009
Teflon Mira Mist, concentric, 35-45 psi, 0.2-2.5ml/min	---	020-060-306
Peek Mira Mist, concentric, 45 psi, 0.2-2.5ml/min	N077-5330	020-060-308

Peristaltic Pump Tubing*

Description	Cross Ref.	Catalog #
Orange/Red, 0.0075", PVC, 2-stops, 12/pak	N069-5476	020-030-002
Black/Black, 0.030", PVC, 2-stops, 12/pak	0990-8587	020-030-007
Red/Red, 0.045", PVC, 2-stops, 12/pak	N068-0375	020-030-010
Black/White, 0.125", PVC, 2-stops, 12/pak	N812-2012	020-030-021
Black/Black, 0.030", PVC, 3-stops, 12/pak	B050-6058	020-033-007
White/White, 0.040", PVC, 3-stops, 12/pak	B050-6742	020-033-009
Red/Red, 0.045", PVC, 3-stops, 12/pak	B019-3160	020-033-010
Purple/Purple, 0.080", PVC, 3-stops, 12/pak	B019-9034	020-033-015
Yellow/Blue/Yellow, 0.060", PVC, 3-stops, 12/pak	B019-3161	020-033-019
Black/White/Black, 0.125", PVC, 3-stops, 12/pak	B050-8310	020-033-021
Black/Black, 0.030", Solvent Flex, 2-stops, 12/pak	0047-3550	020-030-107
Red/Red, 0.045", Solvent Flex, 2-stops, 12/pak	N068-0376	020-030-110
White/White, 0.040", Solvent Flex, 3-stops, 12/pak	B050-7692	020-033-109
Black/Black, 0.030", Silicon, 2-stops, 12/pak	N047-3552	020-030-207
Red/Red, 0.045", Silicon, 2-stops, 12/pak	N069-1595	020-030-210
Black/Black, 0.030", Viton, 2-stops, 12/pak	0990-8587	020-030-307
Perkin Elmer Mixing Kit, PVC, 3-stops, Orange/Orange:Orange/Green, 4/pak	---	020-049-602
Perkin Elmer Sciex Mixing Kit, PVC, 2-stops, Black/Black:Orange/Red, 4/pak	---	020-049-603

* For a complete list, please contact **SCP SCIENCE**

For more information or to place an order, please contact your
SCP SCIENCE Sales Representative or distributor.

USA
Tel.: (800) 361-6820
Fax: (800) 253-5549



Canada / International
Tel.: (800) 361-6820 / (514) 457-0701
Fax: (800) 253-5549 / (514) 457-4499



Europe
Tel.: +33 (0)1 69 18 71 17
Fax: +33 (0)1 60 92 05 67