

SCP SCIENCE Offers a Wide Selection of Standards:

More than 31 elements available at 1000 ppm in 125 and 500 ml bottles for flame and graphite furnace atomic absorption.

- ◆ Certificate of Analysis with actual matrix, concentration and traceability to NIST
- ◆ Complete documentation for audit purposes
- ◆ 2 expiry dates (up to 21 months unopened & 15 months opened)
- ◆ Immediate availability

For Best Results, Use
SCP SCIENCE
AA Single Element
Calibration
Standards

MK-AA-STANDARDS-6.0

ELEMENT	SYMBOL	MATRIX	CODE	CATALOGUE NUMBER	
				125 ml	500 ml
Aluminum	Al	HCl	✓ ⑧	140-002-131	140-002-135
Antimony	Sb	HNO ₃ / tr. Tartaric Acid	✓ ⑧	140-001-511	140-001-515
Arsenic	As	HNO ₃	✓ ⑧	140-001-331	140-001-335
Barium	Ba	HNO ₃	✓ ⑧	140-001-561	140-001-565
Beryllium	Be	HNO ₃	✓ ⑧	140-001-041	140-001-045
Bismuth	Bi	HNO ₃	✓ ⑧	140-001-831	140-001-835
Boron	B	H ₂ O		140-000-051	140-000-055
Cadmium	Cd	HNO ₃	✓ ⑧	140-001-481	140-001-485
Calcium	Ca	HNO ₃	✓ ⑧	140-001-201	140-001-205
Chromium	Cr	HCl	✓ ⑧	140-002-241	140-002-245
Cobalt	Co	HNO ₃	✓ ⑧	140-001-271	140-001-275
Copper	Cu	HNO ₃	✓ ⑧	140-001-291	140-001-295
Gold	Au	HCl	✓ ⑧	140-002-791	140-002-795*
Iron	Fe	HNO ₃	✓ ⑧	140-001-261	140-001-265
Lead	Pb	HNO ₃	✓ ⑧	140-001-821	140-001-825
Lithium	Li	HNO ₃	✓ ⑧	140-001-031	140-001-035
Magnesium	Mg	HNO ₃	✓ ⑧	140-001-121	140-001-125
Manganese	Mn	HNO ₃	✓ ⑧	140-001-251	140-001-255
Mercury	Hg	HNO ₃	✓ ⑧	140-001-801	140-001-805
Molybdenum	Mo	H ₂ O		140-000-421	140-000-425
Nickel	Ni	HNO ₃	✓ ⑧	140-001-281	140-001-285
Potassium	K	HNO ₃	✓ ⑧	140-001-191	140-001-195
Selenium	Se	HNO ₃	✓ ⑧	140-001-341	140-001-345
Silicon	Si	H ₂ O / tr. HF		140-000-141	140-000-145
Silver	Ag	HNO ₃	✓ ⑧	140-001-471	140-001-475
Sodium	Na	HNO ₃	✓ ⑧	140-001-111	140-001-115
Strontium	Sr	HNO ₃	✓ ⑧	140-001-381	140-001-385
Tin	Sn	HCl	✓ ⑧	140-002-501	140-002-505
Titanium	Ti	H ₂ O / tr. HF		140-000-221	140-000-225
Vanadium	V	HNO ₃	✓ ⑧	140-001-231	140-001-235
Zinc	Zn	HNO ₃	✓ ⑧	140-001-301	140-001-305



With every purchase of glassware, receive a FREE **PlasmaCAL** ICP-AES/MS 2 x 25 ml box, 1000µg/ml standards

For a Complete Listing of **SCP SCIENCE** Products, Get your FREE Catalogue Today!!



Catalogue No: 140-000-015



Catalogue No: 140-000-005

Order Online At:
www.scpscience.com
Email:
sales@scpscience.com

North America / International

Tel.: (800) 361-6820
Fax: (800) 253-5549

Europe

Tel.: +33 (0) 1 69 18 71 17
Fax: +33 (0) 1 60 92 05 67

China

Tel.: +86 (10) 58032301
Fax: +86 (10) 58032302

Brazil

Tel.: +55 19 8194-8488