

# SCP SCIENCE

EASY AS

1

2


3

PARTICULATE FREE

## THE NEXT GENERATION IN FILTERING SAMPLES HAS ARRIVED!

The vacuum assisted *Digi*FILTER assembly provides a quick and easy way to filter samples prior to analysis


**4** POSITION  
MANIFOLD



AVAILABLE IN

**100 ml**

FORMAT!



**12** POSITION  
MANIFOLD

**50 ml**



**VIEW VIDEO DEMONSTRATION ONLINE AT**

[www.scpscience.com/en/products/categories?id=382&name=digifilters--accessories](http://www.scpscience.com/en/products/categories?id=382&name=digifilters--accessories)



CORPORATE HQ

CANADA / USA

FRANCE

GERMANY

UK

INDIA

Phone: +1 (514) 457-0701  
Fax: +1 (514) 457-4499

Phone: +1 (800) 361-6820  
Fax: +1 (800) 253-5549


Phone: +33 (0) 1 69 18 71 17  
Fax: +33 (0) 1 60 92 05 67

Phone: +49 (0) 8342-89560-61  
Fax: +49 (0) 8342-89560-69

Phone: +44 (800) 368 6147

Phone: +91 9790975795

**DigiFILTERs** increase lab productivity by filtering particulates from digested samples in a matter of seconds. Available in 0.45, 0.7, 1.0 and 10 µm porosity membranes. Choose the proper membrane type (Teflon®, Metricel, and glass fiber) for your application.



50 ML FORMAT	QTY.	CATALOG #
DigiFILTER, 0.2 µm Hydrophilic Teflon®	100	010-500-088
DigiFILTER, 0.45 µm Hydrophilic Teflon®	100	010-500-070
DigiFILTER, 1.0 µm Hydrophilic Teflon®	100	010-500-080
DigiFILTER, 10.0 µm Hydrophilic Teflon®	100	010-500-074
DigiFILTER, 0.45 µm High Purity Hydrophilic Teflon®	100	010-500-370
DigiFILTER, 1.0 µm High Purity Hydrophilic Teflon®	100	010-500-380
DigiFILTER, 0.45 µm Metricel	100	010-500-072
DigiFILTER, 0.45 µm Mixes Cellulose Estel	100	010-500-118
DigiFILTER, 0.7 µm Glass Fiber for TCLP	100	010-500-082
DigiFILTER, Manifold, 4 position	1	010-500-234
DigiFILTER, Manifold, 12 position	1	010-500-230
DigiTUBEs, RackLock with caps	750	010-500-261
DigiTUBEs, RackLock no caps	750	010-500-262
DigiTUBEs, Non-RackLock with caps	750	010-500-263
DigiTUBEs, Non-RackLock no caps	750	010-500-264

100 ML FORMAT	QTY.	CATALOG #
DigiFILTER, 0.2 µm Hydrophilic Teflon®	50	010-501-088
DigiFILTER, 0.45 µm Hydrophilic Teflon®	50	010-501-070
DigiFILTER, 1.0 µm Hydrophilic Teflon®	50	010-501-080
DigiFILTER, 10.0 µm Hydrophilic Teflon®	50	010-501-074
DigiFILTER, 0.45 µm High Purity Hydrophilic Teflon®	50	010-501-370
DigiFILTER, 1.0 µm High Purity Hydrophilic Teflon®	50	010-501-380
DigiFILTER, 0.7 µm Glass Fiber for TCLP	50	010-501-082
DigiFILTER, Manifold, 4 position	1	010-501-234
DigiFILTER, Manifold, 12 position	1	010-501-238
DigiTUBEs, RackLock with caps	300	010-501-261
DigiTUBEs, RackLock no caps	300	010-501-262
DigiTUBEs, Non-RackLock with caps	300	010-501-263
DigiTUBEs, Non-RackLock no caps	300	010-501-264

VACUUM PUMP	CATALOG #
Standard (115V) 20 L/min.	010-500-235
Standard (115V) 60 L/min.	010-500-237
Standard (230V) 17 L/min.	010-790-509
Standard (230V) 58 L/min.	010-500-239
Teflon® (115V) 20 L/min.	010-790-530
Teflon® (230V) 58 L/min.	010-790-531



Catalog #: 010-500-235 / 237

Catalog #: 010-790-530 / 531

CORPORATE HQ	CANADA / USA	FRANCE	GERMANY	UK	INDIA
Phone: +1 (514) 457-0701 Fax: +1 (514) 457-4499	Phone: +1 (800) 361-6820 Fax: +1 (800) 253-5549	Phone: +33 (0) 1 69 18 71 17 Fax: +33 (0) 1 60 92 05 67	Phone: +49 (0) 8342-89560-61 Fax: +49 (0) 8342-89560-69	Phone: +44 (800) 368 6147	Phone: +91 9790975795