

# DigiFILTER

## 1, 2, 3... Particulate Free!

*The Next Generation in filtering samples has arrived.*

*The vacuum assisted DigiFILTER assembly provides a quick and easy way to filter samples prior to analysis.*



**NEW!**

*\*patent pending*

## SCP SCIENCE

[www.scpscience.com](http://www.scpscience.com)

### USA

Tel.: (800) 361-6820

Fax: (800) 253-5549

### Canada / International

Tel.: (800) 361-6820 /+1 (514) 457-0701

Fax: (800) 253-5549 /+1 (514) 457-4499

### Europe

Tel.: +33 (0) 1 69 18 71 17

Fax: +33 (0) 1 60 92 05 67

### China

Tel.: +86 (10) 58032301


Fax: +86 (10) 58032302

4, 5, 6... Disposable, easy and quick

# DigiFILTER

Filter 12 samplers at a time using the DigiFILTER manifold

Our patent-pending DigiFILTERs will increase lab productivity by filtering particulates from digested samples in a matter of seconds. Available in 0.45 or 1.0 µm hydrophillic Teflon® membranes.



Item	Catalog Number
DigiFILTER, 0.45 micron Teflon® (pk/100)	010-500-070
DigiFILTER, 1.0 micron Teflon® (pk/100)	010-500-080

## DigiFILTER Accessories

Item	Cat. Number
DigiTUBEs, RackLock (750/Case)	010-500-261
DigiTUBEs, Non-RackLock (750/Case)	010-500-263
Vacuum Pump, Standard (115V)	010-500-235
Vacuum Pump, Teflon® lined (115V)	010-790-530
Vacuum Pump, Standard (230V)	010-790-509
Vacuum Pump, Teflon® lined (230V)	010-790-531

Item	Catalog Number
DigiFILTER Manifold, 12 position	010-500-230
DigiFILTER Manifold, 4 position	010-500-234
18 Position Rack	010-500-025



Catalogue Number:  
010-500-070

Catalogue Number:  
010-500-080