

Consistent Even Heating

SCP SCIENCE provides the full range of supplies for all your Graphite Furnace needs. Made from high purity, high-density graphite, these tubes exhibit an extremely low "Coefficient of Thermal Expansion" (CTE), ensuring less stress placed upon the pyrolytic coating, thereby increasing their usable lifetime. Direct replacements for most AA instruments (Perkin Elmer®, Varian®, TJA®, Hitachi® and other OEMs).

- Tubes meet & exceed OEM specifications
- Provides high sensitivity and low noise
- Unconditional money back guaranteed



**SPECIAL
OFFER**



FREE AA STANDARDS*

With the purchase of each pack of graphite furnace tubes receive a 125 ml, 1000 µg/ml AA standard at NO cost!

**Not Included: Supplemental hazardous material charges: Au*

Tube and platform CTE's have been closely matched, reducing the possibility of the platforms releasing from the walls. Thermal conductivity is uniform to allow consistent, even heating, not only within the tube, but from tube-to-tube as well.

www.scpscience.com

USA

Tel.: (800) 361-6820
Fax: (800) 253-5549

Canada / International

Tel.: (800) 361-6820 / +1 (514) 457-0701
Fax: (800) 253-5549 / +1 (514) 457-4499

Europe

Tel.: +33 (0) 1 69 18 71 17
Fax: +33 (0) 1 60 92 05 67

China

Tel.: +86 (10) 58032301
Fax: +86 (10) 58032302

PerkinElmer®			
Description	Quantity (/pk)	Cross Reference	Catalog Number
L'vov Platform Set (pre-inserted)	10	B011-2660 (pk/10), B300-0343 (pk/20)	030-111-001
Platform Tube, Pyrolytically Coated	10	B013-7111 (pk/5), B012-1092 (pk/10), B300-1254 (pk/20), B010-9322 (pk/50)	030-111-002
L'vov Platform	10	B013-7112 (pk/5), B012-1091 (pk/10), B300-1256 (pk/20), B010-9324 (pk/50)	030-111-003
Standard Tube, Pyrolytically Coated	10	B013-5197 (pk/5), B013-5653 (pk/10), B300-0342 (pk/20), B009-1504 (pk/50)	030-111-004
Standard Tube, Uncoated	10	B013-7113 (pk/5), B300-1253 (pk/20), B007-0699 (pk/50)	030-111-005
HGA Contact Set w/Sensor Hole	SET	B012-8490 (pk/1), B018-0363 (pk/5)	030-111-006
Zeeman Contacts	SET	B011-6823 (pk/1), B018-0361 (pk/5)	030-111-007
HGA Contact Set w/o Sensor Hole	SET	B012-8495 (pk/1), B313-0086 (pk/5)	030-111-008
Tube, Pyrolytic w/Fork Ptl.	10	B050-5057 (pk/20)	030-111-010

Shimadzu®			
Description	Quantity (/pk)	Cross Reference	Catalog Number
Standard Tube, Pyrolytically Coated, 60° contact cone	10	200-54525R1	030-117-001
Standard Tube, Uncoated, 60° contact cone	10	200-54520R1	030-117-002
Tube with built-in Platform, Pyrolytically coated, 60° contact cone	10	206-82541R1	030-117-003
Standard Tube, Uncoated, 60° contact cone	10	200-54520R1	030-117-005
Standard Tube, Pyrolytically Coated, 60° contact cone	10	200-54525R1	030-117-004
Tube for Extended Injec. Vol., Uncoated, 60° contact cone	10	206-50587	030-117-007
Tube for Extended Injec. Vol., Pyrolytically Coated, 60° contact cone	10	206-69984-02	030-117-006
Standard Tube, Pyrolytically Coated, 90° contact cone	10	206-50588	030-117-009
Tube with built-in Platform, Pyrolytically coated, 90° contact cone	10	205-50887	030-117-008

Hitachi®			
Description	Quantity (/pk)	Cross Reference	Catalog Number
Standard Tube, Pyrolytically Coated	10	180-7444, 180-7403	030-113-001
Standard Tube, Uncoated	10	180-7400	030-113-002
Forced Platform Tube, Pyrolytically Coated	10	190-6007	030-113-003
Forced Platform	10	190-6008	030-113-006
Standard Platform	10	180-7404	030-113-007
Tube for Extended Injec. Vol., Pyrolytically Coated	10	190-6003	030-113-010
Forced Platform Set (pre-inserted)	10	190-0028	030-113-011
Contacts	1	180-7401	030-113-013

Varian®			
Description	Quantity (/pk)	Cross Reference	Catalog Number
Partition Tube, Pyrolytically Coated	10	63-100012-00	030-110-001
Partition Tube, Pyrolytically Coated, High Purity	10	63-100012-HP	030-110-011
Partition Tube, Pyrolytically Coated, Extended Life	10	63-100012-EL	030-110-021
Platform Tube, NEW!	10	63-100026-00	030-110-025
Platform Tube, Pyrolytically Coated	10	63-100027-00	030-110-002
Bone Platform Tube, Pyrolytically Coated	10	63-100011-00	030-110-003
Forced Platform Tube, Pyrolytically Coated	10	63-100023-00	030-110-012
Forced Platform	10	63-100024-00	030-110-013
Forced Platform Set (pre-inserted)	10	63-100023-90	030-110-004
Monosnap Platform Set (pre-inserted)	2	190032700	030-110-005
Electrode Contact Shroud	1	63-100016-00	030-110-006
Zeeman Electrode Contact	2	63-100018-00	030-110-007
Zeeman Shroud	1	63-100017-00	030-110-008
		63-100019-00	030-110-009

Thermo - Unicam®			
Description	Quantity (/pk)	Cross Reference	Catalog Number
Standard Tube, Uncoated	10	9423 393 90031	030-116-001
Standard Tube, Pyrolytically Coated	10	9423 393 90091	030-116-002
Barrel Tube, Uncoated	10	9423 390 95031	030-116-003
Smooth Tube, Pyrolytically Coated	10	9423 393 95091	030-116-004
Partridge Tube, Extended Life (ELC)	10	9423 393 95041	030-116-005
Partridge Tube, Uncoated	10	9423 393 95031	030-116-006
Partridge Tube, Pyrolytically Coated	10	9423 393 90191	030-116-007
Forced Tube Set (pre-inserted)	10	9423 393 95191	030-116-008
Tube with Silk, Pyrolytically Coated	10	9423 393 95081	030-116-009
Partridge Tube (Smooth), Pyrolytically Coated	10	9423 393 95071	030-116-010
Zeeman Electrode Contact	2	9423 393 95161	030-116-115
Electrode Contact	2	9423 393 95011	030-116-116
Probe	10	9423 393 90081	030-116-119

Thermo - TJA®			
Description	Quantity (/pk)	Cross Reference	Catalog Number
Standard Tube, Uncoated	10	124544-03	030-118-001
Standard Tube, Pyrolytically Coated	10	124544-02	030-118-002
Forced Platform Set (pre-inserted), ADM	10	130941-00	030-118-003
Forced Platform Tube, Pyrolytically Coated	10	---	030-118-004
Graphite Door Kit	1	124778-00	030-118-010
Graphite Cell Body	1	124818-00	030-118-011
Contact Electrodes	2	---	030-118-015

SCP SCIENCE

Providing Innovative Solutions to Analytical Chemists

sales@scpscience.com

Perkin Elmer®, Varian®, Hitachi®, Thermo-TJA®, GBC®, Thermo-Unicam® and Shimadzu® are registered trademarks. SCP SCIENCE is not affiliated or licensed with them.

GBC®			
Description	Quantity (/pk)	Cross Reference	Catalog Number
Standard Tube, Pyrolytically Coated	10	9900 5900	030-112-001
Platform	10	9900 6000	030-112-002
Electrode Contacts	2	9900 6100	030-112-003
GBC Shroud	1	45 0004 00	030-112-005
Platform Set (pre-inserted)	10	---	030-112-006
Ultra Z Tube, Pyrolytically Coated	10	45 0012 00	030-112-007

VOLUME DISCOUNTS AVAILABLE