

XRF Cells



SCP SCIENCE offers Products & Supplies for Liquid, Powder & Solid Sample Analysis. XRF Cells are available in 32, 40, 47 mm. Direct replacement for OEM cells for Oxford, Asoma, Horiba, Philips and ARL.




- Integrated overflow reservoir to collect expanding samples
- Secure snap-on closure preventing spills
- Made from high purity polyethylene
 - Resistant to most acids, bases and most organic chemicals
- Special Discount Offered

XRF Cells

Catalog #	Style	Description	Dia. (mm)	Hght. (mm)	Vol. (mm)	Qty.
040-080-010	Single Open Ended	Snap-Post & Reservoir	32	23	7	200
040-080-011	Single Open Ended	Snap-Post & Reservoir	40	23	13	100
040-080-012	Single Open Ended	Venting	32	23	7	200
040-080-015	Single Open Ended	Trimless	32	25	9	200
040-080-017	Single Open Ended	Trimless	40	25	16	100
040-080-020	Single Open Ended (Horiba)	Trimless	47	19	15	100
040-080-022	Double Open Ended (Horiba)	Trimless	47	19	15	100
040-080-030	Double Open Ended	Trimless	32	25	9	200
040-080-031	Double Open Ended	Trimless	40	25	16	100
040-080-032	Double Open Ended	---	40	35	26	100
040-080-033	Double Open Ended	Provision for Cap	32	23	12	200
040-080-034	Double Open Ended	Cap for 040-080-033	---	---	---	200
040-080-036	Single Open Ended	Reservoir & Venting	40	23	12	100
040-080-038	Single Open Ended	Trimless	45	42	8	100
040-080-040	Double Open Ended	Snap-Post & Reservoir	32	23	7	100
040-080-042	Double Open Ended	Venting	40	23	18	100
040-080-044	Double Open Ended	Snap-Post & Reservoir	40	23	12	100
040-080-046	Oxford Cell System*	Replace aluminum holder	40-	35	26	100
040-080-048	Disposable Cell Liner and caps (Oxford)	---	32	---	---	100
040-080-100	Vent Hole Punch	---	---	---	---	1

* Does not work with autosampler

For more information or to receive your free samples, contact your
SCP SCIENCE Sales Representative today!



Sample
Prep

ICP-AES & MS

AA

XRF

Calibration
Standards

Certified
Reference
Materials

Available!
Oil and Metallo-Organic
Standards



- Full range of concentrations
- Complete with a Certificate of Analysis



MK-MKG003-XRFT-2.0-E

USA

Tel.: (800) 361-6820
Fax: (800) 253-5549

Canada / International

Tel.: (800) 361-6820 / (514) 457-0701
Fax: (800) 253-5549 / (514) 457-4499

Europe

Tel.: +33 (0)1 69 18 71 17
Fax: +33 (0)1 60 92 05 67


XRF Thin Film

Continuous rolls, pre-cut sheets and easy snap sheets are available in different chemical compositions and uniform thickness. Minimum trace levels of impurities ensure reproducible day-to-day results. When selecting a film, resistance to chemical attack, transmission characteristics, and purity must be balanced for optimal results with a specific sample type.

Thin Film - Pre-Cut Sheets

Catalog #	Description	Style	Thickness (mil.)	Dimensions (mm)	Qty.
040-070-055	Microporous	Circles	1.0/25.0	76	100
040-070-105	Mylar®	Circles	0.24/6.0	64	1000
040-070-245	Prolene®	Circles	0.16/4.0	64	1000
040-070-110	Mylar®	Squares	0.14/3.6	64 x 76	500
040-070-120	Mylar®	Squares	0.24/6.0	64 x 76	500
040-070-130	Kapton®	Squares	0.30/7.5	64 x 76	500
040-070-140	Polypropylene	Squares	0.24/6.0	64 x 76	500
040-070-145	Prolene®	Squares	0.16/4.0	64 x 76	500
040-070-147	Prolene®	Easy Snap	0.16/4.0	76	100
040-070-149	Mylar®	Easy Snap	0.14/3.6	76	100
040-070-153	Etnom - S®	Easy Snap	0.12/3.0	64	100
040-070-155	Etnom - S®	Easy Snap	0.12/3.0	76	100

Note: Hostaphan® is an alternate name for Mylar®
Different types of Easy Snap are available on request.



Easy Snap consist of a thin film attached to a carrier frame for easy handling.

Purchase 5 rolls of Mylar® or Kapton® Film and receive 10% off!



Thin Film - Continuous Roll

Catalog #	Description	Thickness (mil./microns)	Width (mm)	Length (m)
040-070-010	Mylar®	0.12/3.0	73	115
040-070-020	Mylar®	0.24/6.0	73	115
040-070-035	Mylar®	0.14/3.6	76	91
040-070-037	Mylar®	0.24/6.0	76	91
040-070-030	Kapton®	0.30/7.5	73	18.8
040-070-040	Polypropylene	0.24/6.0	76	91
040-070-042	Polypropylene	0.50/12.0	76	1
040-070-045	Prolene®	0.16/4.0	76	91
040-070-050	Microporous	1.00/2.5	64	5.1
040-070-100	Film Dispenser*	---	---	---

* For Mylar and Kapton film only

For more information, contact your
SCP SCIENCE Sales Representative today!