

SCP SCIENCE PVC 2-stop regular Pump Tubing Cross-Reference Guide

	ID (mm)	JY - HORIBA	PerkinElmer	SPECTRO	Teledyne	Thermo - iCAP, Elemental, TJA	Thermo - VG	Varian (Agilent)
022-030-002	0.19		N0695476					
022-030-003	0.25		N0773117					
022-030-004	0.38		N0777110	47201037				
022-030-005	0.51				309000691	1161920		
022-030-006	0.64		N8122012	47101001		842312051511		
022-030-007	0.76	46.791.022	09908587		3093552			3710027200
022-030-008	0.89			47101003	3093551			
022-030-009	1.02	46.791.022				842312051521	1200212	3710034400
022-030-010	1.14		N0680375		30900694			
022-030-011	1.30	46.791.018		47501010				3710034500
022-030-012	1.42		N0777521					
022-030-013	1.65			47101014				3710034600
022-030-014	1.85	46.791.016	N0777518					
022-030-015	2.06		N0777561					
022-030-016	2.29							3710027300
022-030-017	2.54							
022-030-018	2.79							
022-030-019	1.52		09923536	47101006				
022-030-021	3.18		N8122012					
022-030-022	1.75			47101005				
022-030-023	0.95							G182065220
022-030-024	0.44		N0773113					