

FOR IMMEDIATE RELEASE

New Product from **SCP SCIENCE** – *NOVAWAVE* Microwave Digestion System

(Montreal, February 17th, 2010) - **SCP SCIENCE** is pleased to announce the introduction of the *NOVAWAVE* microwave digestion system.

NOVAWAVE delivers a dramatic improvement in sample throughput and speed of digestion plus a major reduction in labour and human intervention compared to any other digestion system. A microwave sample digestion tunnel system employing 12 dynamically operating micro cavities to process 12 samples simultaneously. Other features include individual sample temperature control and monitoring; seamless microwave power control with no duty cycle; built-in EPA and DIN methods; and USB/Ethernet port connectivity for data processing and storage.



Up to 168 samples at once – With the *NOVAWAVE* Transporter attachment, prepare up to 14 racks of 12 samples at a time. The Transporter acts as a conveyor, automatically loading a new rack of samples into the microwave and removing the digested samples to the auto-cooling and auto-venting stations.

Versatility – *NOVAWAVE* is available in two configurations: stand alone Model SA and fully automated Model FA; the former can be upgraded at any time.

About SCP SCIENCE

Founded in 1980, **SCP SCIENCE** is a successful privately owned manufacturer and distributor of analytical equipment, supplies, standards, reagents, and certified reference materials for the inorganic analytical laboratories market. The company has multiple manufacturing units for instruments, glass and quartz ware, injection moulding and much more plastic ware, aqueous and oil based standards. The headquarters of **SCP SCIENCE** are located in Montreal, Quebec, Canada. The company maintains additional facilities in Canada, the USA, and Europe.

Contact:

Carl Chan - Marketing Manager
SCP SCIENCE – www.scpscience.com
Tel: (514) 457-0701 • Fax: (514) 457-4499 • cchan@scpscience.com

U.S.A

348 Route 11
Champlain NY
12919-4816
Tel: 1-800-361-6820
Fax: 1-800-253-5549

CANADA

21800 Clark Graham
Baie Durfé, QC
H9X 4B6
Tel: +1 (514) 457-0701
Fax: +1 (514) 457-4499

EUROPE

12 avenue du Québec
Bâtiment 1-2, Silic 642
91965 Courtaboeuf CIDEX
Tel: +33 (0) 1 69 18 71 17
Fax: +33 (0) 1 60 92 05 67

CHINA

No 2. Building Room 727
106 Ma Jia Bao Road
Feng Tai District. Beijing 100068
Tel.: +86 (10) 58032301
Fax: +86 (10) 58032302