

 ORE RESEARCH & EXPLORATION P/L
 ABN 28 006 859 856

 37A Hosie Street · Bayswater North · VIC 3153 · AUSTRALIA

 ▲ 61 3 9729 0333

 ⑤ info@ore.com.au

 ● www.ore.com.au

CERTIFICATE OF ANALYSIS FOR

LATERITIC NICKEL-COBALT ORE CERTIFIED REFERENCE MATERIAL OREAS 181



Table 1. Certified Values	, SDs, 9	95% Confidence and	Tolerance Limits	for OREAS 181.
---------------------------	----------	--------------------	-------------------------	----------------

Constituent	Certified	100	95% Confid	dence Limits	95% Tolerance Limits					
Constituent	Value	150	Low	High	Low	High				
Borate Fusion XRF										
Al ₂ O ₃ , Aluminium(III) oxide (wt.%)	11.56	0.118	11.49	11.62	11.50	11.61				
CaO, Calcium oxide (wt.%)	2.37	0.020	2.36	2.38	2.35	2.39				
Co, Cobalt (ppm)	451	11	446	456	435	466				
Cr ₂ O ₃ , Chromium(III) oxide (wt.%)	1.24	0.016	1.23	1.25	1.23	1.25				
Cu, Copper (ppm)	74	13	66	82	IND	IND				
Fe ₂ O ₃ , Iron(III) oxide (wt.%)	35.94	0.205	35.85	36.04	35.79	36.09				
K ₂ O, Potassium oxide (wt.%)	0.132	0.005	0.129	0.135	0.130	0.135				
MgO, Magnesium oxide (wt.%)	2.05	0.041	2.03	2.08	2.03	2.08				
MnO, Manganese oxide (wt.%)	0.168	0.005	0.166	0.171	0.166	0.171				
Na ₂ O, Sodium oxide (wt.%)	0.448	0.020	0.436	0.459	0.431	0.465				
Ni, Nickel (ppm)	5123	59	5091	5155	5071	5175				
P ₂ O ₅ , Phosphorus(V) oxide (wt.%)	0.017	0.003	0.015	0.018	IND	IND				
SiO ₂ , Silicon dioxide (wt.%)	33.79	0.114	33.74	33.83	33.60	33.97				
SO ₃ , Sulphur trioxide (wt.%)	0.189	0.008	0.181	0.197	0.182	0.196				
TiO ₂ , Titanium dioxide (wt.%)	0.339	0.008	0.335	0.343	0.331	0.347				
Zn, Zinc (ppm)	96	7.8	91	102	IND	IND				
Thermogravimetry										
LOI ¹⁰⁰⁰ , Loss On Ignition @1000°C	11.26	0.339	11.06	11.46	11.20	11.32				
(wt.%)										
	- 44.04	0.050	11.10	11.50	11.00	44.50				
Al2O3, Aldminium(III) Oxide (wt. %)	11.31	0.352	11.10	11.52	11.09	11.52				
Ba, Banum (ppm)	135	9	130	139	129	141				
	2.37	0.106	2.31	2.43	2.32	2.42				
	< 10					10.0				
Ce, Cenum (ppm)	12.3	1.2	11.0	13.7	11.9	12.8				
Co, Cobait (ppm)	451	10	447	454	441	461				
Cr ₂ O ₃ , Chromium(III) oxide (wt.%)	1.23	0.043	1.20	1.25	1.20	1.25				
Cs, Cesium (ppm)	0.57	0.035	0.55	0.59	0.52	0.62				
Cu, Copper (ppm)	//	4.4	/4	80	73	81				
	0.43	80.0	0.34	0.52						
Fe ₂ O ₃ , Iron(III) oxide (wt.%)	35.34	0.895	34.83	35.86	34.71	35.97				
Ga, Gallium (ppm)	10.5	0.79	9.7	11.2	9.7	11.2				
Gd, Gadolinium (ppm)	1.35	0.21	1.12	1.57	1.20	1.49				
Hf, Hafnium (ppm)	1.74	0.22	1.46	2.02	IND	IND				
Ho, Holmium (ppm)	0.27	0.03	0.24	0.29	IND	IND				
K ₂ O, Potassium oxide (wt.%)	0.125	0.015	0.116	0.133	0.109	0.141				
La, Lanthanum (ppm)	4.04	0.58	3.42	4.66	3.73	4.35				
Li, Lithium (ppm)	12.3	1.6	10.7	13.9	IND	IND				
MgO, Magnesium oxide (wt.%)	2.01	0.062	1.97	2.05	1.97	2.05				
MnO, Manganese oxide (wt.%)	0.166	0.007	0.162	0.169	0.160	0.172				
Mo, Molybdenum (ppm)	< 5	IND	IND	IND	IND	IND				
Na ₂ O, Sodium oxide (wt.%)	0.420	0.036	0.380	0.460	0.404	0.436				

Note: intervals may appear asymmetric due to rounding



	Certified		95% Confid	dence Limits	95% Tolerance Limits							
Constituent	Value	1SD	Low	High	Low	High						
Borate / Peroxide Fusion ICP continued												
Nb, Niobium (ppm)	1.96	0.38	1.63	2.30	IND	IND						
Nd, Neodymium (ppm)	5.30	0.53	4.87	5.72	4.68	5.91						
Ni, Nickel (ppm)	5048	126	5012	5084	4914	5182						
Pr, Praseodymium (ppm)	1.29	0.17	1.15	1.43	IND	IND						
Rb, Rubidium (ppm)	5.61	0.384	5.38	5.84	5.35	5.87						
S, Sulphur (wt.%)	0.071	0.006	0.067	0.074	IND	IND						
Sb, Antimony (ppm)	1.19	0.16	1.10	1.28	IND	IND						
Sc, Scandium (ppm)	38.2	1.27	37.4	38.9	36.5	39.8						
SiO ₂ , Silicon dioxide (wt.%)	33.32	1.050	32.64	34.01	32.63	34.02						
Sm, Samarium (ppm)	1.36	0.22	1.15	1.57	IND	IND						
Sr, Strontium (ppm)	87	5.1	84	90	85	89						
Tb, Terbium (ppm)	0.22	0.04	0.20	0.25	IND	IND						
Th, Thorium (ppm)	2.72	0.128	2.60	2.83	2.58	2.85						
TiO ₂ , Titanium dioxide (wt.%)	0.333	0.015	0.325	0.341	0.322	0.344						
U, Uranium (ppm)	0.98	0.073	0.93	1.03	0.82	1.14						
V, Vanadium (ppm)	146	15	135	157	137	154						
W, Tungsten (ppm)	3.22	0.38	2.95	3.50	IND	IND						
Y, Yttrium (ppm)	5.73	0.64	5.39	6.07	5.47	6.00						
Yb, Ytterbium (ppm)	0.89	0.080	0.84	0.94	IND	IND						
Zn, Zinc (ppm)	100	10	91	108	92	107						
Zr, Zirconium (ppm)	61	4.9	57	66	59	64						
Infrared Combustion	<u>.</u>											
C, Carbon (wt.%)	0.653	0.028	0.635	0.671	0.640	0.665						
S, Sulphur (wt.%)	0.059	0.012	0.052	0.067	0.055	0.064						
Gas / Liquid Pycnometry												
SG, Specific Gravity (Unity)	2.87	0.155	2.74	3.00	2.82	2.91						

Table 1 continued.

Note: intervals may appear asymmetric due to rounding

INTRODUCTION

OREAS reference materials are intended to provide a low cost method of evaluating and improving the quality of analysis of geological samples. To the geologist they provide a means of implementing quality control in analytical data sets generated in exploration from the grass roots level through to prospect evaluation, and in grade control at mining operations. To the analyst they provide an effective means of calibrating analytical equipment, assessing new techniques and routinely monitoring in-house procedures.

SOURCE MATERIALS

OREAS 181 has been prepared from a blend of low and high grade Ni-Co lateritic ores sourced from the Bulong deposit located 35km east of Kalgoorlie in Western Australia. This ore formed from prolonged lateritic weathering of Archaean (Yilgarn Craton) olivine rich ultramafic/komatiite flows. Grades of >1% Ni were generated in zones of more intense



weathering associated with faulting and bedrock alteration. The Ni-Co nontronitic (Fe-Ni smectite clays) siliceous ores at Bulong formed with a goethitic overprint in the upper laterite profile.

COMMINUTION AND HOMOGENISATION PROCEDURES

The material constituting OREAS 181 was prepared in the following manner:

- drying to constant mass at 105°C;
- crushing and milling to 100% minus 35 microns;
- preliminary homogenisation and check assaying of low and high grade source materials;
- final homogenisation by blending the source materials in specific ratios to achieve target grades;
- packaging in 10g units sealed in laminated foil pouches and 1kg units in plastic jars.

ANALYTICAL PROGRAM

Seventeen commercial analytical laboratories participated in the program to certify the analytes reported in Table 1. The following methods were employed:

- Borate fusion with XRF for common nickel laterite assemblage (up to 14 laboratories depening on the analyte);
- Thermogravimetric analysis of LOI at 1000°C (15 labs);
- Borate or peroxide fusion for full elemental suite ICP-OES and ICP-MS finishes (up to 15 laboratories depending on the element);
- C and S by IR combustion furnace (14 labs);
- Specific gravity by gas (6 labs) or liquid (3 labs) pycnometry.

For the round robin program ten 250g test units were taken at predetermined intervals during the bagging stage, immediately following homogenisation and are considered representative of the entire 270kg batch. The six samples received by each laboratory were obtained by taking one 20g scoop split from each of six different test units. This format maximised representivity of the parent batch at each lab. Table 1 presents the 63 certified values together with their associated 1SD's, 95% confidence and tolerance limits and Table 2 below shows 26 indicative values. Table 3 provides performance gate intervals for the certified values based on their associated pooled standard deviations. Tabulated results of all elements together with analytical method codes, uncorrected means, medians, standard deviations, relative standard deviations and per cent deviation of lab means from the corrected mean of means (PDM³) are presented in the detailed certification data for this CRM (**OREAS 181 DataPack.xIsx**).



							1		
Constituent	Unit	Value	Constituent	Unit	Value	Constituent	Unit	Value	
Borate Fusion XRF									
As	ppm	27.8	Pb	ppm	< 50	Zr	ppm	78	
BaO	ppm	157	Sc	ppm	38.7				
CI	ppm	699	V_2O_5	ppm	247				
Borate / Peroxide	Fusion ICP								
Ag	ppm	< 5	Ge	ppm	1.85	Sn	ppm	1.67	
As	ppm	23.9	In	ppm	< 0.1	Та	ppm	0.16	
В	ppm	75	Lu	ppm	0.14	Те	ppm	< 1	
Be	ppm	< 1	P_2O_5	wt.%	0.019	TI	ppm	< 0.5	
Bi	ppm	0.075	Pb	ppm	9.23	Tm	ppm	0.13	
Dy	ppm	1.22	Re	ppm	< 0.1				
Er	ppm	0.78	Se	ppm	< 20				

Table 2. Indicative Values for OREAS 181.

Note: the number of significant figures reported is not a reflection of the level of certainty of stated values. They are instead an artefact of ORE's in-house CRM-specific LIMS.

STATISTICAL ANALYSIS

Certified Values, Confidence Limits, Standard Deviations and Tolerance Limits (Table 1) have been determined for each analyte following removal of individual, laboratory dataset (batch) and 3SD outliers (single iteration). For individual outliers within a laboratory batch the z-score test is used in combination with a second method that determines the per cent deviation of the individual value from the batch median. Outliers in general are selected on the basis of z-scores > 2.5 and with per cent deviations (i) > 3 and (ii) more than three times the average absolute per cent deviation for the batch. In certain instances statistician's prerogative has been employed in discriminating outliers. Each laboratory data set mean is tested for outlying status based on z-score discrimination and rejected if > 2.5. After individual and laboratory data set (batch) outliers have been eliminated a non-iterative 3 standard deviation filter is applied, with those values lying outside this window also relegated to outlying status. The Certified Values are the means of accepted laboratory means after outlier filtering.

The 95% Confidence Limits are inversely proportional to the number of participating laboratories and inter-laboratory agreement. It is a measure of the reliability of the certified value. A 95% confidence interval indicates a 95% probability that the true value of the analyte under consideration lies between the upper and lower limits. *95% Confidence Limits should not be used as control limits for laboratory performance.*

Standard Deviation values (1SDs) are reported in Table 1 and provide an indication of a level of performance that might reasonably be expected from a laboratory being monitored by this CRM in a QA/QC program. The SD's take into account errors attributable to measurement uncertainty and CRM variability. For an effective CRM the contribution of the latter should be negligible in comparison to measurement errors. The SD values thus include all sources of measurement uncertainty: between-lab variance, within-run variance (precision errors) and CRM variability. OREAS prepared reference materials have a level of homogeneity such that the observed variance from repeated analysis has its origin almost exclusively in the analytical process rather than the reference material itself.



The SD for each analyte's certified value is calculated from the same filtered data set used to determine the certified value, i.e. after removal of any individual, lab dataset (batch) and 3SD outliers (single iteration). These outliers can only be removed after the absolute homogeneity of the CRM has been independently established, i.e. the outliers must be confidently deemed to be analytical rather than arising from inhomogeneity of the CRM. The standard deviation is then calculated for each analyte from the pooled accepted analyses generated from the certification program.

In the application of SD's in monitoring performance it is important to note that not all laboratories function at the same level of proficiency and that different methods in use at a particular laboratory have differing levels of precision. Each laboratory has its own inherent SD (for a specific concentration level and analyte-method pair) based on the analytical process and this SD is not directly related to the round robin program.

The majority of data generated in the round robin program was produced by a selection of world class laboratories. The SD's thus generated are more constrained than those that would be produced across a randomly selected group of laboratories. To produce more generally achievable SD's the 'pooled' SD's provided in this report include inter-lab bias. This 'one size fits all' approach may require revision at the discretion of the QC manager concerned following careful scrutiny of QC control charts.

Table 3 shows **Performance Gates** calculated for two and three standard deviations. As a guide these intervals may be regarded as warning or rejection for multiple 2SD outliers, or rejection for individual 3SD outliers in QC monitoring, although their precise application should be at the discretion of the QC manager concerned. A second method utilises a 5% window calculated directly from the certified value. Standard deviation is also shown in relative percent for one, two and three relative standard deviations (1RSD, 2RSD and 3RSD) to facilitate an appreciation of the magnitude of these numbers and a comparison with the 5% window. Caution should be exercised when concentration levels approach lower limits of detection of the analytical methods employed as performance gates calculated from standard deviations tend to be excessively wide whereas those determined by the 5% method are too narrow.

Tolerance Limits (ISO Guide 3207) were determined using an analysis of precision errors method and are considered a conservative estimate of true homogeneity. The meaning of tolerance limits may be illustrated for nickel (Ni) by fusion XRF, where 99% of the time (1- α =0.99) at least 95% of subsamples (p=0.95) will have concentrations lying between 5071 and 5175 ppm. Put more precisely, this means that if the same number of subsamples were taken and analysed in the same manner repeatedly, 99% of the tolerance intervals so constructed would cover at least 95% of the total population, and 1% of the tolerance intervals would cover less than 95% of the total population (ISO Guide 35). *Please note that tolerance limits pertain to the homogeneity of the CRM only and should not be used as control limits for laboratory performance.*

Based on the statistical analysis of the results of the inter-laboratory certification program it can be concluded that OREAS 181 is fit-for-purpose as a certified reference material (see 'Intended Use' below).



Constituent	Constituent Certified Value	Absolute Standard Deviations				Relative Standard Deviations			5% window		
Constituent		1SD	2SD Low	2SD High	3SD Low	3SD High	1RSD	2RSD	3RSD	Low	High
Borate Fusion	Borate Fusion XRF										
Al ₂ O ₃ , wt.%	11.56	0.118	11.32	11.79	11.20	11.91	1.02%	2.04%	3.06%	10.98	12.13
CaO, wt.%	2.37	0.020	2.33	2.41	2.31	2.43	0.84%	1.68%	2.51%	2.25	2.49
Co, ppm	451	11	429	472	418	483	2.41%	4.81%	7.22%	428	473
Cr ₂ O ₃ , wt.%	1.24	0.016	1.21	1.27	1.19	1.29	1.32%	2.64%	3.95%	1.18	1.30
Cu, ppm	74	13	47	100	34	113	17.89%	35.79%	53.68%	70	77
Fe ₂ O ₃ , wt.%	35.94	0.205	35.53	36.35	35.33	36.56	0.57%	1.14%	1.71%	34.15	37.74
K ₂ O, wt.%	0.132	0.005	0.122	0.143	0.116	0.148	4.00%	8.00%	12.01%	0.126	0.139
MgO, wt.%	2.05	0.041	1.97	2.14	1.93	2.18	1.98%	3.97%	5.95%	1.95	2.16
MnO, wt.%	0.168	0.005	0.158	0.179	0.153	0.184	3.02%	6.04%	9.07%	0.160	0.177
Na ₂ O, wt.%	0.448	0.020	0.407	0.488	0.387	0.509	4.53%	9.06%	13.59%	0.425	0.470
Ni, ppm	5123	59	5004	5242	4945	5301	1.16%	2.32%	3.48%	4867	5379
P ₂ O ₅ , wt.%	0.017	0.003	0.011	0.022	0.008	0.025	16.77%	33.53%	50.30%	0.016	0.018
SiO ₂ , wt.%	33.79	0.114	33.56	34.01	33.44	34.13	0.34%	0.68%	1.02%	32.10	35.47
SO ₃ , wt.%	0.189	0.008	0.174	0.204	0.166	0.212	4.04%	8.08%	12.12%	0.179	0.198
TiO ₂ , wt.%	0.339	0.008	0.323	0.354	0.315	0.362	2.31%	4.62%	6.93%	0.322	0.356
Zn, ppm	96	7.8	81	112	73	120	8.07%	16.15%	24.22%	92	101
Thermogravin	netry										
LOI ¹⁰⁰⁰ , wt.%	11.26	0.339	10.59	11.94	10.25	12.28	3.01%	6.02%	9.02%	10.70	11.83
Borate / Pero	xide Fusion	ICP									
Al ₂ O ₃ , wt.%	11.31	0.352	10.60	12.01	10.25	12.37	3.11%	6.23%	9.34%	10.74	11.87
Ba, ppm	135	9	117	152	108	161	6.52%	13.05%	19.57%	128	141
CaO, wt.%	2.37	0.106	2.16	2.58	2.05	2.69	4.47%	8.94%	13.41%	2.25	2.49
Cd, ppm	< 10	IND	IND	IND	IND	IND	IND	IND	IND	IND	IND
Ce, ppm	12.3	1.2	9.9	14.8	8.6	16.1	10.10%	20.20%	30.30%	11.7	13.0
Co, ppm	451	10	430	472	420	482	2.32%	4.63%	6.95%	428	473
Cr ₂ O ₃ , wt.%	1.23	0.043	1.14	1.31	1.10	1.35	3.47%	6.94%	10.42%	1.16	1.29
Cs, ppm	0.57	0.035	0.50	0.64	0.46	0.67	6.20%	12.40%	18.60%	0.54	0.60
Cu, ppm	77	4.4	68	86	64	90	5.66%	11.32%	16.98%	73	81
Eu, ppm	0.43	0.08	0.27	0.59	0.19	0.67	18.80%	37.61%	56.41%	0.41	0.45
Fe ₂ O ₃ , wt.%	35.34	0.895	33.55	37.13	32.66	38.03	2.53%	5.06%	7.59%	33.57	37.11
Ga, ppm	10.5	0.79	8.9	12.0	8.1	12.8	7.58%	15.16%	22.74%	9.9	11.0
Gd, ppm	1.35	0.21	0.93	1.76	0.72	1.97	15.54%	31.07%	46.61%	1.28	1.41
Hf, ppm	1.74	0.22	1.30	2.18	1.07	2.41	12.77%	25.54%	38.31%	1.65	1.83
Ho, ppm	0.27	0.03	0.21	0.33	0.18	0.35	10.92%	21.85%	32.77%	0.25	0.28
K ₂ O, wt.%	0.125	0.015	0.095	0.155	0.080	0.170	12.04%	24.08%	36.12%	0.119	0.131
La, ppm	4.04	0.58	2.88	5.20	2.30	5.78	14.36%	28.71%	43.07%	3.84	4.24
Li, ppm	12.3	1.6	9.2	15.5	7.6	17.0	12.80%	25.60%	38.39%	11.7	12.9
MgO, wt.%	2.01	0.062	1.89	2.13	1.82	2.19	3.07%	6.14%	9.21%	1.91	2.11
MnO, wt.%	0.166	0.007	0.153	0.179	0.146	0.186	3.99%	7.98%	11.98%	0.158	0.174
Mo, ppm	< 5	IND	IND	IND	IND	IND	IND	IND	IND	IND	IND
Na ₂ O, wt.%	0.420	0.036	0.348	0.492	0.312	0.528	8.59%	17.19%	25.78%	0.399	0.441
Nb, ppm	1.96	0.38	1.19	2.73	0.81	3.11	19.59%	39.18%	58.76%	1.86	2.06
Nd, ppm	5.30	0.53	4.23	6.36	3.70	6.89	10.06%	20.12%	30.18%	5.03	5.56

Table 3. Performance Gates for OREAS 181.

Note: intervals may appear asymmetric due to rounding.



Constituent	Certified	Absolute Standard Deviations				Relative Standard Deviations			5% window		
Constituent	Value	1SD	2SD Low	2SD High	3SD Low	3SD High	1RSD	2RSD	3RSD	Low	High
Borate / Peroxide Fusion ICP continued											
Ni, ppm	5048	126	4796	5300	4670	5426	2.50%	5.00%	7.50%	4796	5300
Pr, ppm	1.29	0.17	0.96	1.62	0.80	1.79	12.77%	25.54%	38.31%	1.23	1.36
Rb, ppm	5.61	0.384	4.84	6.38	4.46	6.76	6.84%	13.68%	20.51%	5.33	5.89
S, wt.%	0.071	0.006	0.059	0.082	0.053	0.088	8.42%	16.85%	25.27%	0.067	0.074
Sb, ppm	1.19	0.16	0.87	1.50	0.72	1.66	13.24%	26.49%	39.73%	1.13	1.25
Sc, ppm	38.2	1.27	35.6	40.7	34.4	42.0	3.32%	6.64%	9.96%	36.3	40.1
SiO ₂ , wt.%	33.32	1.050	31.22	35.42	30.17	36.47	3.15%	6.30%	9.45%	31.66	34.99
Sm, ppm	1.36	0.22	0.93	1.79	0.71	2.01	15.92%	31.85%	47.77%	1.29	1.43
Sr, ppm	87	5.1	77	98	72	103	5.87%	11.74%	17.61%	83	92
Tb, ppm	0.22	0.04	0.14	0.31	0.09	0.36	19.67%	39.34%	59.01%	0.21	0.24
Th, ppm	2.72	0.128	2.46	2.97	2.33	3.10	4.71%	9.42%	14.13%	2.58	2.85
TiO ₂ , wt.%	0.333	0.015	0.304	0.363	0.289	0.378	4.44%	8.89%	13.33%	0.317	0.350
U, ppm	0.98	0.073	0.83	1.13	0.76	1.20	7.48%	14.95%	22.43%	0.93	1.03
V, ppm	146	15	116	175	102	190	10.12%	20.24%	30.36%	139	153
W, ppm	3.22	0.38	2.46	3.99	2.08	4.37	11.83%	23.67%	35.50%	3.06	3.38
Y, ppm	5.73	0.64	4.44	7.02	3.80	7.67	11.24%	22.48%	33.72%	5.45	6.02
Yb, ppm	0.89	0.080	0.73	1.05	0.65	1.13	8.98%	17.96%	26.94%	0.84	0.93
Zn, ppm	100	10	79	120	69	130	10.18%	20.36%	30.54%	95	105
Zr, ppm	61	4.9	51	71	46	76	8.04%	16.08%	24.11%	58	64
Infrared Com	bustion										
C, wt.%	0.653	0.028	0.597	0.708	0.570	0.736	4.25%	8.50%	12.76%	0.620	0.685
S, wt.%	0.059	0.012	0.035	0.084	0.022	0.097	20.94%	41.89%	62.83%	0.057	0.062
Gas / Liquid F	Pycnometry										
SG, Unity	2.87	0.155	2.56	3.18	2.40	3.33	5.40%	10.80%	16.20%	2.73	3.01

Table 3 continued.

Note: intervals may appear asymmetric due to rounding.

PARTICIPATING LABORATORIES

- 1. ALS, Brisbane, QLD, Australia
- 2. ALS, Lima, Peru
- 3. ALS, Vancouver, BC, Canada
- 4. Argile Analytica, Calgary, Alberta, Canada
- 5. Bureau Veritas Commodities Canada Ltd, Vancouver, BC, Canada
- 6. Bureau Veritas Geoanalytical, Adelaide, SA, Australia
- 7. Bureau Veritas Geoanalytical, Perth, WA, Australia
- 8. Inspectorate (BV), Lima, Peru
- 9. Intertek Genalysis, Perth, WA, Australia
- 10. Intertek Testing Services, Cupang, Muntinlupa, Philippines
- 11. Intertek Testing Services, Townsville, QLD, Australia
- 12. Nagrom, Perth, WA, Australia
- 13. Ni Lab, Pouembout, New Caledonia
- 14. PT Geoservices Ltd, Cikarang, Jakarta Raya, Indonesia
- 15. SGS Australia Mineral Services, Perth, WA, Australia
- 16. SGS Geosol Laboratorios Ltda, Vespasiano, Minas Gerais, Brazil
- 17. SGS Mineral Services, Townsville, QLD, Australia



PREPARER AND SUPPLIER

Certified reference material OREAS 181 is prepared, certified and supplied by:



ORE Research & Exploration Pty LtdTel:+613-9729 033337A Hosie StreetFax:+613-9729 8338Bayswater North VIC 3153Web:www.ore.com.auAUSTRALIAEmail:info@ore.com.au

It is available in unit sizes of 10g (single-use laminated foil pouches) and 1kg (plastic jars).

INTENDED USE

OREAS 181 is intended for the following uses:

- for the monitoring of laboratory performance in the analysis of analytes reported in Table 1 in geological samples;
- for the verification of analytical methods for analytes reported in Table 1;
- for the calibration of instruments used in the determination of the concentration of analytes reported in Table 1.

STABILITY AND STORAGE INSTRUCTIONS

OREAS 181 is an oxidised reference material and is stable in the laminated foil pouches. Under normal conditions of storage it has a shelf life beyond ten years.

INSTRUCTIONS FOR CORRECT USE

The certified values determined via fusion ICP, C and S by infrared combustion furnace and SG by pycnometry refer to the concentration levels in the packaged state. There is no need for drying prior to weighing and analysis.

In contrast the certified values determined via borate fusion XRF and for LOI at 1000° C are on a dry basis. This requires the removal of hygroscopic moisture by drying in air to constant mass at 105° C. If the reference material is not dried prior to analysis, the certified values should be corrected to the moisture-bearing basis.

HANDLING INSTRUCTIONS

Fine powders pose a risk to eyes and lungs and therefore standard precautions such as the use of safety glasses and dust masks are advised.

TRACEABILITY

The analytical samples were selected in a manner to represent the entire batch of prepared CRM. This 'representivity' was maintained in each submitted laboratory sample



batch and ensures the user that the data is traceable from sample selection through to the analytical results that underlie the consensus values. Each analytical data set has been validated by its assayer through the inclusion of internal reference materials and QC checks during analysis. The laboratories were chosen on the basis of their competence (from past performance in inter-laboratory programs) for a particular analytical method, analyte or analyte suite, and sample matrix. Most of these laboratories have and maintain ISO 17025 accreditation. The certified values presented in this report are calculated from the means of accepted data following robust statistical treatment as detailed in this report.

LEGAL NOTICE

Ore Research & Exploration Pty Ltd has prepared and statistically evaluated the property values of this reference material to the best of its ability. The Purchaser by receipt hereof releases and indemnifies Ore Research & Exploration Pty Ltd from and against all liability and costs arising from the use of this material and information.

QMS ACCREDITED

ORE Pty Ltd is accredited to ISO 9001:2008 by Lloyd's Register Quality Assurance Ltd for its quality management system including development, manufacturing, certification and supply of CRMs.



CERTIFYING OFFICER

Craig Hamlyn (B.Sc. Hons - Geology), Technical Manager - ORE P/L

REFERENCES

ISO Guide 30 (1992), Terms and definitions used in connection with reference materials.

ISO Guide 31 (2000), Reference materials – Contents of certificates and labels.

ISO Guide 3207 (1975), Statistical interpretation of data - Determination of a statistical tolerance interval.

ISO Guide 35 (2006), Certification of reference materials - General and statistical principals.

