

SCP SCIENCE

New colored caps for *DigiTUBEs*

**Differentiate samples
with different colored caps**

Red



Yellow



Orange



Clear



USA

Tel.: (800) 361-6820
Fax: (800) 253-5549

Canada / International

Tel.: (800) 361-6820 /+1 (514) 457-0701
Fax: (800) 253-5549 /+1 (514) 457-4499

Europe

Tel.: +33 (0) 1 69 18 71 17
Fax: +33 (0) 1 60 92 05 67

China

Tel.: +86 (10) 51905671
Fax: +86 (10) 51905672

New Caps



NEW

Catalogue No.	Description	Quantity
010-500-060	Blue caps for <i>DigiTUBE</i> s	(pk/250)
010-500-130	Yellow caps for <i>DigiTUBE</i> s	(pk/250)
010-500-140	Orange caps for <i>DigiTUBE</i> s	(pk/250)
010-500-150	Red caps for <i>DigiTUBE</i> s	(pk/250)
010-500-172	Clear caps for <i>DigiTUBE</i> s	(pk/250)

DigiTUBE containers replace and eliminate the need for five different laboratory vessels (i.e digestion beaker, graduated cylinder, volumetric flask, autosampler tube, and storage container). Manufactured using virgin polypropylene to avoid metal contamination. For use in AA or ICP analysis sample storage, and other applications.



RackLOCK style

RackLOCK feature allows one-handed closure of screw cap digestion tubes saving time.

Material: Polypropylene (PP)

- Dimensions: 30mm x 110mm (62ml)
- Graduated: 5, 10, 15, 20, 25, 30, 35, 40, 45, 50 ml
Class A tolerance at the 50 ml line

Ordering Information

Catalogue No.	Description	Quantity
010-500-261	<i>DigiTUBE</i> - RackLock, with blue caps	750
010-500-262	<i>DigiTUBE</i> - RackLock, without caps, choose your colored caps	750
010-500-263	<i>DigiTUBE</i> - Non-RackLock, with blue caps	750
010-500-264	<i>DigiTUBE</i> - Non-RackLock, without caps, choose your colored caps	750