

XRF Cells



SCP SCIENCE offers Products & Supplies for Liquid, Powder & Solid Sample Analysis. XRF Cells are available in 32, 40, 47 mm. Direct replacement for OEM cells for Oxford, Asoma, Horiba, Philips and ARL.

Ask about our volume discounts of 20% - 30% off!




- Integrated overflow reservoir to collect expanding samples
- Secure snap-on closure preventing spills
- Made from high purity polyethylene
 - Resistant to most acids, bases and most organic chemicals
- Special Discount Offered

XRF Cells

Catalog #	Style	Description	Dia. (mm)	Hght. (mm)	Vol. (mm)	Qty.
040-080-010	Single Open Ended	Snap-Post & Reservoir	32	23	7	200
040-080-011	Single Open Ended	Snap-Post & Reservoir	40	23	13	100
040-080-012	Single Open Ended	Venting	32	23	7	200
040-080-015	Single Open Ended	Trimless	32	25	9	200
040-080-017	Single Open Ended	Trimless	40	25	16	100
040-080-020	Single Open Ended (Horiba)	Trimless	47	19	15	100
040-080-022	Double Open Ended (Horiba)	Trimless	47	19	15	100
040-080-030	Double Open Ended	Trimless	32	25	9	200
040-080-031	Double Open Ended	Trimless	40	25	16	100
040-080-032	Double Open Ended	---	40	35	26	100
040-080-033	Double Open Ended	Provision for Cap	32	23	12	200
040-080-034	Double Open Ended	Cap for 040-080-033	---	---	---	200
040-080-036	Single Open Ended	Reservoir & Venting	40	23	12	100
040-080-038	Single Open Ended	Trimless	45	42	8	100
040-080-040	Double Open Ended	Snap-Post & Reservoir	32	23	7	100
040-080-042	Double Open Ended	Venting	40	23	18	100
040-080-044	Double Open Ended	Snap-Post & Reservoir	40	23	12	100
040-080-046	Oxford Cell System*	Replace aluminum holder	40-	35	26	100
040-080-048	Disposable Cell Liner and caps (Oxford)	---	32	---	---	100
040-080-100	Vent Hole Punch	---	---	---	---	1

* Does not work with autosampler

For more information or to receive your free samples, contact your SCP SCIENCE Sales Representative today!



www.scpscience.com

Sample Prep

ICP-AES & MS

AA

XRF

Calibration Standards

Certified Reference Materials

Available!
Oil and Metallo-Organic Standards



- Full range of concentrations
- Complete with a Certificate of Analysis

Call For More Information!

USA

Tel.: (800) 361-6820
Fax: (800) 253-5549

Canada / International

Tel.: (800) 361-6820 / (514) 457-0701
Fax: (800) 253-5549 / (514) 457-4499

Europe

Tel.: +33 (0)1 69 18 71 17
Fax: +33 (0)1 60 92 05 67


XRF Thin Film

Continuous rolls, pre-cut sheets and easy snap sheets are available in different chemical compositions and uniform thickness. Minimum trace levels of impurities ensure reproducible day-to-day results. When selecting a film, resistance to chemical attack, transmission characteristics, and purity must be balanced for optimal results with a specific sample type.

Thin Film - Pre-Cut Sheets

Catalog #	Description	Style	Thickness (mil.)	Dimensions (mm)	Qty.
040-070-055	Microporous	Circles	1.0/25.0	76	100
040-070-105	Mylar®	Circles	0.24/6.0	64	1000
040-070-245	Prolene®	Circles	0.16/4.0	64	1000
040-070-110	Mylar®	Squares	0.14/3.6	64 x 76	500
040-070-120	Mylar®	Squares	0.24/6.0	64 x 76	500
040-070-130	Kapton®	Squares	0.30/7.5	64 x 76	500
040-070-140	Polypropylene	Squares	0.24/6.0	64 x 76	500
040-070-145	Prolene®	Squares	0.16/4.0	64 x 76	500
040-070-147	Prolene®	Easy Snap	0.16/4.0	76	100
040-070-149	Mylar®	Easy Snap	0.14/3.6	76	100
040-070-153	Etnom - S®	Easy Snap	0.12/3.0	64	100
040-070-155	Etnom - S®	Easy Snap	0.12/3.0	76	100

Note: Hostaphan® is an alternate name for Mylar®
Different types of Easy Snap are available on request.



Easy Snap consist of a thin film attached to a carrier frame for easy handling.

Purchase 5 rolls of Mylar® or Kapton® Film and receive 10% off!



Thin Film - Continuous Roll

Catalog #	Description	Thickness (mil./microns)	Width (mm)	Length (m)
040-070-010	Mylar®	0.12/3.0	73	115
040-070-020	Mylar®	0.24/6.0	73	115
040-070-035	Mylar®	0.14/3.6	76	91
040-070-037	Mylar®	0.24/6.0	76	91
040-070-030	Kapton®	0.30/7.5	73	18.8
040-070-040	Polypropylene	0.24/6.0	76	919
040-070-042	Polypropylene	0.50/12.0	76	1
040-070-045	Prolene®	0.16/4.0	76	91
040-070-050	Microporous	1.00/2.5	64	5.1
040-070-100	Film Dispenser*	---	---	---

* For Mylar and Kapton film only

For more information, contact your
SCP SCIENCE Sales Representative today!