


## AGILENT USERS

OPTIMIZE YOUR SAMPLE INTRODUCTION WITH

# SCP SCIENCE DEMOUNTABLE AND SEMI-DEMOUNTABLE TORCHES FOR ICP SPECTROMETRY

- ▶ The semi-demountable and demountable torches, with removable injectors, allow faster transition between different sample matrices
- ▶ Variety of quartz injector sizes (0.8 mm, 1.4 mm, 1.8 mm, 2.4 mm) with replaceable extended outer tube, with slot
- ▶ Dual View (DV) or Radial View (RV) according to your needs
- ▶ The Easy-Fit design automatically aligns the torch and connects gases for speed and convenience



PART NUMBER	DESCRIPTION *	REPLACES AGILENT	ALSO REPLACES AGILENT
020-056-160	5100 Series, Fully Demountable, with 0.8 mm Quartz Injector, RV	G8020-68001	
020-056-158	5100 Series, Fully Demountable with 1.4 mm Quartz Injector, RV	G8020-68007	
020-056-156	5100 Series, Fully Demountable, with 1.8 mm Alumina Injector, RV		G8010-60234
020-056-159	5100 Series, Fully Demountable with 1.4 mm Quartz Injector, DV	G8020-68002	
020-056-155	5100 Series, Fully Demountable, with 1.8 mm Quartz Injector, DV	G8020-68005	
020-056-157	5100 Series, Fully Demountable with 1.8 mm Alumina Injector, DV	G8020-68003	G8010-60231
020-056-163	5100 Series, Fully Demountable with 2.4 mm Quartz Injector, DV	G8020-68004	G8010-60232
020-056-151	5100 Series, Semi-Demountable, with 1.4 mm Quartz Injector, DV	G8010-60237	G8012-60000
020-056-153	5100 Series, Semi-Demountable, with 1.8 mm Quartz Injector, DV	G8010-60236	G8010-60228
020-056-060	5100 Series, 0.8 mm Quartz Injector		
020-056-058	5100 Series, 1.4 mm Quartz Injector		
020-056-054	5100 Series, 1.8 mm Quartz Injector		
020-056-062	5100 Series, 2.4 mm Quartz Injector		
020-056-056	5100 Series, 1.8 mm Alumina Injector		
020-056-020	5100 Series, Quartz Outer Tube set, DV		
020-056-022	5100 Series, Quartz Outer Tube set, RV		

\*Replacement injectors for use with SCP SCIENCE Agilent torches only (catalogue no. 020-056-1xx)

PROMO

**FREE** pack of *PlasmaFLOW* PVC pump tubing or **FREE** ICP Standard\*

(Single element, 1,000 ppm, 2 x 25 ml) with the purchase of this *PlasmaFLOW* torch.

\*The following elements are excluded from all special offers: Be, Sc, Ru, Rh, Pd, Re, Os, Ir, Pt, Au, Tm, Lu.



CORPORATE HQ	CANADA / USA	FRANCE	GERMANY	UK	INDIA
Phone: +1 (514) 457-0701 Fax: +1 (514) 457-4499	Phone: +1 (800) 361-6820 Fax: +1 (800) 253-5549	Phone: +33 (0) 1 69 18 71 17 Fax: +33 (0) 1 60 92 05 67	Phone: +49 (0) 8342-89560-61 Fax: +49 (0) 8342-89560-69	Phone: +44 (800) 368 6147	Phone: +91 9790975795