

Autosampler Tubes and Caps

12 ml and 15 ml

SCP Science autosampler tubes are manufactured from virgin polypropylene (PP) with ultra-low metal content, making them ideal for applications involving ICP or ICP-MS analysis. These round bottom tubes are available in two convenient sizes: 16 x 100 mm (12 ml) and 17 x 100 mm (15 ml) and designed to meet the requirements of position racks used by most instrument manufacturers.



Key Features

- ▶ Stable in the presence of common acids, alkalis and solvents
- ▶ Made with easy-to-read, Class A graduation lines with ± 0.10 ml tolerances
- ▶ 8 different colors caps available to help identify samples

Meet the demands of your daily lab routines with **SCP SCIENCE** autosampler tubes.

Description	Quantity	Part Number
16 x 100 mm (12 ml) Polypropylene autosampler tubes	1000 tubes	010-515-609
	2 bags of 1000 tubes	010-515-629
17 x 100 mm (15 ml) Polypropylene autosampler tubes	1000 tubes	010-515-607
	2 bags of 1000 tubes	010-515-627

Color caps manufactured from Low-Density Polyethylene (LDPE) with ultra-low metal content. All caps are made with a double seal flange to guarantee a leak-free seal with the side wall of the tube.

Description	Quantity	Part Number
White	1000 caps	010-515-620
Yellow	1000 caps	010-515-630
Orange	1000 caps	010-515-640
Red	1000 caps	010-515-650
Blue	1000 caps	010-515-660
Clear	1000 caps	010-515-672
Green	1000 caps	010-515-680
Purple	1000 caps	010-515-690