

SCP SCIENCE introduces a bench-top fume hood for *DigiPREP MS, Jr.*, and *HP* sample digestion systems. Replacement air is drawn through the dual HEPA filter inlets offering cleaner ambient air. A powerful blow unit constantly pulls air into the fume hood and exhausts it out of the building.



MK-FUMEHOOD39-2-0-E

- ◆ Integrated Lexan door for easy access of samples
- ◆ Filtered air with HEPA filter adaptors for cleaner air flow
- ◆ Manufactured with Lexan offering durable and see-through shell constructions: Assembled in minutes and requires no tools

Ordering Information & Specifications

DigiPREP Fume Hood System

Catalog Number	010-500-039
Construction	Corian bottom, Lexan walls
Size (in)	26W x 23.5D x 20H
Weight (lb)	49

Options & Accessories

DESCRIPTION	QUANTITY	CATALOGUE NUMBER
Variable speed blower unit (115V/230V)	1	010-500-033
• 3 foot exhaust hose with 90 degree elbow*	1	
• 6 foot exhaust hose**	1	
10" Exhaust Hose Extension	1	010-500-034
HEPA Filter Adaptors***	2	010-500-050
Universal Scrubber System	1	010-520-160

* To connect hood and blower unit

** To connect blower unit to fumehood or fume duct

*** Removes 99.97 percent of .3 micron particles and greater from the air

Use the Universal Scrubber System with your DigiPREP Fume Hood!



The Universal Scrubber System removes up to 98% of all generated acid fumes:

- Robust vacuum pump system
- All acid resistant tubing
- Free up fumehood space and eliminate using "water-jet" vacuum system

Call for more information or to request a quotation!!

CANADA • USA

Tel.: (800) 361-6820 / +1 (514) 457-0701
Fax: (800) 253-5549 / +1 (514) 457-4499

FRANCE

Tel.: +33 (0) 1 69 18 71 17
Fax: +33 (0) 1 60 92 05 67

INTERNATIONAL

Tel.: +1 (514) 457-0701
Fax: +1 (514) 457-4499

Order Online At:

www.scpscience.com
Email: sales@scpscience.com