


Whether it's for ICP-AES, ICP-MS or AAS, **SCP SCIENCE** has the autosampler cups and test tubes you need!!

Our Sample Cups and Tubes Are:

- ◆ Competitively priced.
- ◆ An exact match (equivalent) to OEM Products: 100% Satisfaction Guaranteed.
- ◆ Molded either from premium grade polystyrene for clarity, or polypropylene for impact strength and stability to common acids, alkalis, and solvents.

Put them on the same shipment as your Standards & Glassware order !!



Ordering Information Call 1-800-361-6820 to place your order (Toll FREE in North America)

Description & Volume	Quantity/ per case	Composition	Equivalent to	Catalogue No.
1.5 - 2.0, Conical	1000	Polystyrene	B0119079	080-070-102
3.0 - 4.0, Conical	1000	Polystyrene	B0129303	080-070-103
3.0 - 4.0	1000	Polyethylene	9910028400	080-070-105
6.0	1000	Polypropylene	B0193235	130-012-001
16.0	1000	Polypropylene	N7301205	130-010-015
17.0	1000	Polypropylene	n/a	130-012-007
50 Centrifuge Tubes with caps	500	Polypropylene	n/a	130-010-051
Plug Stopper for 130-012-015	1000	n/a	n/a	130-010-115
Plug Stopper, 16 mm, for 130-012-016	1000	n/a	n/a	130-012-016
Plug Stopper, 17 mm, for 130-012-007	1000	n/a	n/a	130-012-017
<i>DigiTUBEs</i> with Screw Caps* (50 ml)	750	Polypropylene	n/a	010-500-263

SCP SCIENCE is the best source for ICP and Atomic Absorption Spectroscopy Supplies: Quartzware, Standards & Reagents, Pump Tubing, Hollow Cathode Lamps, Furnace Tubes and more... Order Online Today!

Online: www.scpscience.com
Email: sales@scpscience.com

Have You Used the NEW *DigiTUBE*?

Perfect for Use in Many Autosamplers!!



Cleanest Tubes Available!!

The *DigiTUBE* is Recognized for

- Leak-free screw caps
- The 50 ml graduation with class "A" accuracy
- The racklock system, for ease of use

Now *DigiTUBE* Just Got Better!!

Made from NEW and improved polypropylene, the *DigiTUBE* now has the lowest catalytic/additive metal content in the market today !!

Call for more information on our new volume packaging and our new low volume pricing!!

MK-MKG003-AUTO-3.0-E

USA
 Tel.: (800) 361-6820
 Fax: (800) 253-5549

Canada / International
 Tel.: (800) 361-6820 • +1 (514) 457-0701
 Fax: (800) 253-5549 • +1 (514) 457-4499

Europe
 Tel.: +33 (0) 1 69 18 71 17
 Fax: +33 (0) 1 60 92 05 67

China
 Tel.: +86 (10) 58032301
 Fax: +86 (10) 58032302