


DigiVAC *Bringing your digestions out of the Fume Hood*

SCP SCIENCE presents an alternative to digesting samples in the conventional fume hood. Use the *DigiVAC* Accessory Set with a *DigiPREP* to move the digestion out of the fume hood and free up valuable space.



Benefits of the DigiVAC


- "Free-up" valuable fume hood space
- Have uniform evaporation rates across all samples
- Observe sample digestions easily through clear hood

Ordering Information DigiVAC system

Part Number	Description
010-500-031	<ul style="list-style-type: none"> • Exhaust hood and corian ring • 2 / 24 position DigiVAC Racks to be used in conjunction
010-500-033	<ul style="list-style-type: none"> • Variable speed blower unit (115V / 230V) • 3 foot exhaust hose with 90 degree elbow to connect hood and blower unit. • 6 foot exhaust hose to connect blower unit to fume hood or fume duct.

DigiVAC options and accessories

010-500-034	10' Exhaust Hose Extension
010-500-027	Extra 24-Position DigiTUBE Racks (DigiVAC version)



Sample
Prep

ICP-AES & MS

AA

XRF

Calibration
Standards

Certified
Reference
Materials

Have you used the
NEW DigiTUBE?



**cleanest tubes
available !!**

DigiTUBE

is recognized for

- leak-free screw caps
- the 50 ml graduation with class "A" accuracy
- the racklock system, for ease of use

Now DigiTUBE just got better !!

- made from NEW and improved polypropylene, the **DigiTUBE** now has the lowest catalytic/additive metal content in the market today !!

Call for more information on our new volume packaging and our new low volume pricing !!

MK-MKGOO3-DVAC-1.0-E

USA

Tel.: (800) 361-6820
Fax: (800) 253-5549

Canada / International

Tel.: (800) 361-6820 / (514) 457-0701
Fax: (800) 253-5549 / (514) 457-4499

Europe

Tel.: +33 (0)1 69 18 71 17
Fax: +33 (0)1 60 92 05 67