

SCP SCIENCE PVC 3-stop regular Pump Tubing Cross-Reference Guide

	ID (mm)	JY - HORIBA	PerkinElmer	SPECTRO	Teledyne	Thermo - iCAP, Elemental, TJA	Thermo - VG	Varian (Agilent)
022-033-002	0.19		N0695477					G183365572
022-033-003	0.25							
022-033-004	0.38					13017706		
022-033-005	0.51							
022-033-006	0.64							
022-033-007	0.76		B0506058				1301670	
022-033-008	0.89					13017700		
022-033-009	1.02		B0506742					
022-033-010	1.14		B0193160			13017701		
022-033-011	1.30							
022-033-012	1.42							
022-033-013	1.65							
022-033-014	1.85							
022-033-015	2.06		B0199034					
022-033-016	2.29							
022-033-017	2.54							
022-033-018	2.79							
022-033-019	1.52		B0193161					
022-033-021	3.18		B0508310					
022-033-022	1.75							
022-033-023	0.95							
022-033-024	0.44							